


ASBA Seat and Captain Seat Frame Mounting Hardware For Elevating Seat



ASBA Seat and Captain Seat Frame Mounting Hardware For Elevating Seat

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1043627-NLA	NLA - Screw, Hex Head w/ Patch (5/16-18 x 1-3/4")	1	4
2		1106521-NLA	NLA - Bracket, Seat Frame Support, Black	1	2
3		7930337-NLA	NLA - Package, Screw, Hex Head w/ Patch (5/16-18 x 1-1/4")	25	4
4		1106522-NLA	NLA - Bracket, Seat Frame Mounting, Black	1	1
5		1039860	Screw, Hex Head w/ Patch (5/16-18 x 3/4")	1	4
6		1026158	Locknut (5/16-18)	1	12
7		1096795-NLA	NLA - Plug Button, Square, Black (1")	1	2
8	A	1108224-NLA	NLA - Hardware, Seat Frame Mounting	1	1

NOTE: A - Includes items 1-7