




Front Riggings - Myon Junior Footrest Support



Front Riggings - Myon Junior Footrest Support

Front Riggings - Myon Junior Hanger

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	8881428752	Hanger Assembly w/ Hardware, 80 Degree, Short - Right	1	1
2			Hanger Bar, 80 Degree, Short, Color - Right	1	1
4			Screw (M6 x 45mm)	1	1
5			Screw (M5 x .8 x 25mm)	1	1
6			Pivot, Hanger, 80 Degree	1	1
7			Cap, Cylinder, Hanger	1	1
8			Screw (M5 x .8 x 23mm)	1	3
9			Screw (M5 x .8 x 40mm)	1	1
10			Collar, Anti-Rattle	1	1
11			Washer (M6 x 12mm x 0.8mm)	1	1
12			Spring, Double	1	1
13			Lock, Latch Release, Female	1	1
14			Lock, Latch Release, Male	1	1
15			Housing, Latch	1	1
16			Bushing, Latch Housing	1	1
17			Nut, Lock (M5 x .8)	1	5
18			Bushing, Latch	1	1
7		0021-20036	Cap, Cylinder, Hanger	1	1
10		1026124	Collar, Anti-Rattle	1	1
12		0572-06603	Spring, Double	1	1
1	A	8881428753	Hanger Assembly w/ Hardware, 80 Degree, Short - Left	1	1
2			Hanger Bar, 80 Degree, Short, Color - Left	1	1
4			Screw (M6 x 45mm)	1	1
5			Screw (M5 x .8 x 25mm)	1	1
6			Pivot, Hanger, 80 Degree	1	1
7			Cap, Cylinder, Hanger	1	1
8			Screw (M5 x .8 x 23mm)	1	3
9			Screw (M5 x .8 x 40mm)	1	1
10			Collar, Anti-Rattle	1	1
11			Washer (M6 x 12mm x 0.8mm)	1	1
12			Spring, Double	1	1
13			Lock, Latch Release, Female	1	1
14			Lock, Latch Release, Male	1	1
15			Housing, Latch	1	1
16			Bushing, Latch Housing	1	1
17			Nut, Lock (M5 x .8)	1	5
18			Bushing, Latch	1	1
10		1026124	Collar, Anti-Rattle	1	1
12		0572-06603	Spring, Double	1	1
1	B	8881514392	Hanger Assembly w/ Hardware, 80 Degree, Long - Right	1	1
3			Hanger Bar, 80 Degree, Long, Color - Right	1	1
4			Screw (M6 x 45mm)	1	1
5			Screw (M5 x .8 x 25mm)	1	1

Front Riggings - Myon Junior Footrest Support

Item Number	Note	Part Number	Description	Quantity Per		
				Package	Assembly	
6	B	0572-06603 8881514393	Pivot, Hanger, 80 Degree	1	1	
7			Cap, Cylinder, Hanger	1	1	
8			Screw (M5 x .8 x 23mm)	1	3	
9			Screw (M5 x .8 x 40mm)	1	1	
10			Collar, Anti-Rattle	1	1	
11			Washer (M6 x 12mm x 0.8mm)	1	1	
12			Spring, Double	1	1	
13			Lock, Latch Release, Female	1	1	
14			Lock, Latch Release, Male	1	1	
15			Housing, Latch	1	1	
16			Bushing, Latch Housing	1	1	
17			Nut, Lock (M5 x .8)	1	5	
18			Bushing, Latch	1	1	
12			Spring, Double	1	1	
1			Hanger Assembly w/ Hardware, 80 Degree, Long - Left	1	1	
3			Hanger Bar, 80 Degree, Long, Color - Left	1	1	
4			Screw (M6 x 45mm)	1	1	
5			Screw (M5 x .8 x 25mm)	1	1	
6		Pivot, Hanger, 80 Degree	1	1		
7		Cap, Cylinder, Hanger	1	1		
8		Screw (M5 x .8 x 23mm)	1	3		
9		Screw (M5 x .8 x 40mm)	1	1		
10		Collar, Anti-Rattle	1	1		
11		Washer (M6 x 12mm x 0.8mm)	1	1		
12		Spring, Double	1	1		
13		Lock, Latch Release, Female	1	1		
14		Lock, Latch Release, Male	1	1		
15		Housing, Latch	1	1		
16		Bushing, Latch Housing	1	1		
17		Nut, Lock (M5 x .8)	1	5		
18		Bushing, Latch	1	1		
10		1026124	Collar, Anti-Rattle	1	1	
12		0572-06603	Spring, Double	1	1	
20		C	1169143	Kit, Hardware, Block Sector	1	1
21				Screw (M6 x 1 x 60mm)	1	1
22				Position Latch	1	1
NOTE: A - Includes items 2, and 4-19. B - Includes items 3, and 4-19. C - Includes items 21 and 22.						