


Cyclical Lever Drive - 7 Spoke



Cyclical Lever Drive - 7 Spoke

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	A,C	1088186-NLA	NLA - Hardware	1	1
1		1025716	Screw, Phillips Pan Head (#10-32 x 11/16")	1	3
2		1087504-NLA	NLA - Plate, Support Rear Outer Hub	1	3
3		1087407-NLA	NLA - Screw, Hex Head (7/16-20 x 6-1/2")	1	1
4		1087354-NLA	NLA - Spacer, Hub Inner (7/16 x 3/4 x 1-1/16")	1	1
5		1084822-NLA	NLA - Washer, Axle (7/16 x 2-3/8" x 1/16")	1	1
6		361441	Cap, End, Notched, Black	1	1
7		1045608	Locknut (7/16-20)	1	1
8		1087908-NLA	NLA - Spacer, Hub Outer Rear (7/16 x 3/4 x 19/32")	1	1
9		1088008	Transmission and Wheel Lock Assembly - Right	1	1
9		1088009	Transmission and Wheel Lock Assembly - Left	1	1
10	D	1087991-NLA	NLA - Steering Link Assembly - Right	1	1
10		1087992-NLA	NLA - Steering Link Assembly - Left	1	1
11		1088144	Fork, Right (6")	1	1
11		1088145	Fork, Left (6")	1	1
11		1088146	Fork, Right (8")	1	1
11		1088147	Fork, Left (8")	1	1
12		B	1088167-NLA	NLA - Kit, CLD 7 Spoke (6" Caster) - Right	1
12	B	1088168-NLA	NLA - Kit, CLD 7 Spoke (6" Caster) - Left	1	1
12	B	1088169-NLA	NLA - Kit, CLD 7 Spoke (8" Caster) - Right	1	1
12	B	1088170-NLA	NLA - Kit, CLD 7 Spoke (8" Caster) - Left	1	1
NOTE: A - Includes items 1 - 8. B - Includes items 1 - 11. C - Not shown. D - Order Item 12 (1088167, 1088168, 1088169, and 1088170) as a replacement for this item.					
Installation of the CLD MUST be performed by a qualified technician. If a non-certified individual performs any work on these units, the Manufacturer's warranty is void.					