


# Shroud - Top

IRC5LXAW and IRC5LXO2AW



## Shroud - Top

### IRC5LXAW and IRC5LXO2AW

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,C	1108803	Cabinet Cover without Tie Wraps	1	1
2		1070050	Tie Wrap (3-29/32" L)	1	1
3		60125500	Cord, Elastic	1	1
4		1115868	Kit, Cabinet Handle	1	1
			Handle, Top, Warm Grey	1	1
5	D	1172003	NLA - Screw, Hex Head Tap (1/4-15 x 3/4")	1	2
6			Package, Power Cord Strap w/ Hardware	1	1
7			Strap, Hook and Loop	10	1
			Rivet, Plastic	10	1
			NLA - Retainer, Push On	10	1
8		1107400	Grid, Filter	1	2
9		1107412	Filter, Cabinet	1	2
10		1110766	Door, Filter Access with Foam Gasket	1	1
	D	1112296-NLA	NLA - Kit, Packaging w/ Carton, Inserts and Owner's Manual	1	1
	B,D	1118389-NLA	NLA - Owner's Manual	1	1

NOTE: A - Includes Cabinet and Items 1-9.  
B - To be used with the IRC5LXAW and IRC5LXO2AW.  
C - After 6/28/07, the Hook and Loop Strap are no longer used on the Cabinet Cover Assembly. If Hook and Loop Strap is needed you must order Package, Power Cord Strap w/ Hardware, Part Number 1172003.  
D - Not Shown.