




Cabinet and Base

Platinum 5 Models



Cabinet and Base Platinum 5 Models

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,E	1108803	Cabinet Cover without Tie Wraps	1	1
2		1006525	Tie Wrap (3-29/32" L)	1	1
3		60125500	Cord, Elastic	1	1
	F	1115868	Kit, Cabinet Handle	1	1
4			Handle, Top, Warm Grey	1	1
5			NLA - Screw, Hex Head Tap (1/4-15 x 3/4")	1	2
6		1107412	Filter, Cabinet	1	1
7		1107400	Grid, Filter	1	1
8		1107414	Foam, Side Shroud	1	2
9		1107413	Foam, Front Shroud	1	1
10	E,J	1172003	Package, Power Cord Strap w/ Hardware	1	1
11			Strap, Hook and Loop	10	1
12			Rivet, Plastic	10	1
			NLA - Retainer, Push On	10	1
13		1110766	Door, Filter Access with Foam Gasket	1	1
14	B	1112259-NLA	NLA - Base Cabinet Assembly	1	1
15		1108816	Foam, Base	1	1
16		1104151	Fitting, Hose Barb with Nut (1/8")	1	1
17		1108211	Cap, End, Dust, Green	1	1
18		1105515-NLA	NLA - Tech Only - Module, Inlet (15A/ 250 VAC)	1	1
19	D,G	1076830-NLA	NLA - Caster (2")	1	4
19		1110426	Caster Assembly (2")	1	4
20		2000576	Stand-Off, Twist-Lok (1/2")	1	1
21		2000360	Grommet, Plastic	1	4
22		2000508	Screw, Phillips Truss Head Self Tap (#8-32 x 5/8")	1	4
24		1062794-NLA	NLA - Cord, Power - Argentina, Australia, Peoples Republic of China and New Zealand	1	1
24	C,F	1062791-NLA	NLA - Cord, Line - Brazil, Korea and Turkey	1	1
24		1131308-NLA	NLA - Tech Only - Cord, Power, Continental Europe	1	1
24		1102634	Cord, Power - India	1	1
24		1071398	Cord, Power - Israel	1	1
24		1062793-NLA	NLA - Cord, Power - Hong Kong, Singapore and Vietnam	1	1
	G,H,J	1071428	Label, Spanish - IRC5LX	1	1
	G,H,J	1071430	Label, Spanish - IRC5LXO2	1	1
	G,H,J	1079322	Label, Brazilian Portuguese - IRC5LX	1	1
	G,H,J	1079323	Label, Brazilian Portuguese - IRC5LXO2	1	1

NOTE: A - Includes Cover Cabinet and Items 2-13.
 B - Includes Base Cabinet and Items 15-21.
 C - Item no longer available, refer to Part #1131308.
 D - Item no longer available, refer to Part #1110426.
 E - After 6/28/07, the Hook and Loop Strap are no longer used on the Cabinet Cover Assembly. If Hook and Loop Strap is needed you must order Package, Power Cord Strap w/ Hardware, Part Number 1172003.
 F - Includes items 4 and 5.
 G - Not shown.
 H - Labels are not available for the following countries: Australia, Peoples Republic of China, Hong Kong, India, Israel, Korea, New Zealand, Singapore, and Vietnam.
 J - Standard equipment manufactured with full English label applied. Language specific overlays are used to customize equipment.