

# Quick Reference Guide - Electric Bed & Pendants

## Pendants

### Full Electric



IVC Pendant

### Semi-Electric




IVC Pendant




### CareGuard



### Pendant Clip



026975



5880M063

# Quick Reference Guide - Electric Bed & Pendants

## Adapter Cords

Bed Type	Bed Model Number	Motor Part Number	Adapter Cord Part Number
Semi-Electric	5310, 5886, and 5887	5880M067-NLA 5880M068-NLA	1104342-NLA See Note(s) A
Semi-Electric	5108-7 5208-7 5308-7	5880M067-NLA 5880M068-NLA	52000M600-NLA See Note(s) B
Full Electric	5410, 5996, and 5997	5880M067-NLA 5990M018-NLA 5880M068-NLA	1104342-NLA See Note(s) A
Full Electric	5409-7 5409-8	5880M067-NLA 5990M018-NLA 5880M068-NLA	52000M600-NLA See Note(s) B

NOTE: A - Connects the White, 4-pin square connector of the new style motor to the white, 6-pin rectangular receiver of the old-style junction box.  
B - Connects the older motors to the new junction box.

### Semi Electric - U.S.

Bed Models	Foot Spring	Plug Type Motor	Pendant (Rectangular Plug)	Junction / Controller Box Kit	Head Motor with Pull Tube	Foot Motor with Pull Tube
5108-7 5208-7 5308-7	5886-NLA	Rectangular	7001M298	52000M368 See Note(s) A	5880M067-NLA See Note(s) B	5880M068-NLA See Note(s) B
5108-8 5208-8 5308-8	5887-NLA	Rectangular	7001M298 4-Way	52000M368 See Note(s) A	5880M067-NLA See Note(s) B	5880M068-NLA See Note(s) B
5110 5210 5310	5890	Square 4 Pin White	7001M298	1104345-NLA	7001M300-NLA	7001M305-NLA
5310IVC	5890IVC	Molded Grey Plug with O-Ring 6 Pin White	1115288	1117177	1121270	1116644

NOTE: A - Junction Box Replacement Kit includes: Junction Box, Power Cord, and Adapter Cords. Adaptor Cords allow the use of rectangular plug motors with square plug Junction Box Kit.  
B - Order corresponding new style motor 5310 Bed and also (1) Adaptor Cord.

### Semi Electric - Canadian

Bed Model	Foot Spring	Plug Type (Motor)	Pendant (Rectangular Plug)	Junction / Controller Box Kit	Head Motor with Pull Tube	Foot Motor with Pull Tube
5311	5891	Square 4 Pin White	7001M298	1104345-NLA	7001M327-NLA	7001M325-NLA
5311IVC	5891IVC	Molded Grey Plug with O-Ring 6 Pin White	1119170-NLA	1119166	1121270	1116644

# **Quick Reference Guide - Electric Bed & Pendants**

---

# Quick Reference Guide - Electric Bed & Pendants

## Semi Electric - International

Bed Models	Foot Spring	Plug Type (Motor)	Pendant (Rectangular Plug)	Junction / Controller Box Kit	Head Motor with Pull Tube	Foot Motor with Pull Tube
5312	5892-NLA	Square 4 Pin White	7001M298	1104346-NLA	7001M324-NLA	7001M325-NLA
ES1000LT	7050LT-NLA	Square 4 Pin White	7001M298	1104345-NLA	7001M303-NLA	7001M307-NLA

## Full Electric and BAR600 - U.S.

Bed Model	Foot Spring	Plug Type (Motor)	Pendant (Rectangular Plug)	Controller/Junction Box Kit	Head Motor with Pull Tube	Hi/Lo Motor with Bracket	Foot Motor with Pull Tube
5409-7	5996-NLA	Rectangular	7001M297	52000M367 See Note(s) A	5880M067-NLA See Note(s) B	5990M018-NLA See Note(s) B	5880M068-NLA See Note(s) B
5409-8	5997-NLA	Rectangular	7001M297	52000M367 See Note(s) A	5880M067-NLA See Note(s) B	5990M018-NLA See Note(s) B	5880M068-NLA See Note(s) B
5410	5490-NLA	Square 4 Pin White	7001M297	1104343	7001M300-NLA	7001M308-NLA	7001M305-NLA
5410IVC	5490IVC	Molded Grey Plug with O-Ring 6 Pin White	1115290	1117178	1121270	1116653	1116644
BAR600	BAR5490	Square 4 Pin White	7001M297	1104347-NLA See Note(s) C	1125307	1132750	5600M008
BAR600IVC	BAR5490IVC	Square 4 Pin White	7001M297	1106398-NLA See Note(s) C	1125307	1132750	5600M008

NOTE: A - Junction Box Replacement Kit includes Junction Box, Power Cord, and Adapter Cords. Adapter Cords allow the use of Rectangular Plug Motors with current Square Plug Junction Box.  
 B - Order corresponding new style Motor from 5410 Bed and Adapter Cord (1) ea. p/n 1104342.  
 C - No longer available, use Kit part #1104343 used for model 5410 as its replacement. Includes the Controller/Junction Box and Installation Instructions.

## Full Electric - Canadian

Bed Model	Foot Spring	Plug Type (Motor)	Pendant (Rectangular Plug)	Junction / Controller Box Kit	Head Motor with Pull Tube	Hi/Lo Motor with Bracket	Foot Motor with Pull Tube
5411	5491	Square 4 Pin White	7001M297	1104343	7001M327-NLA	7001M329	7001M328-NLA

# Quick Reference Guide - Electric Bed & Pendants

## Full Electric - Canadian

Bed Model	Foot Spring	Plug Type (Motor)	Pendant (Rectangular Plug)	Junction / Controller Box Kit	Head Motor with Pull Tube	Hi/Lo Motor with Bracket	Foot Motor with Pull Tube
5411IVC	5491IVC	Molded Grey Plug with O-Ring 6 Pin White	1119171	1119167	1121270	1116653	1116644

# Quick Reference Guide - Electric Bed & Pendants

## Full Electric - International

Bed Model	Foot Spring	Plug Type (Motor)	Pendant (Rectangular Plug)	Junction / Controller Box Kit	Head Motor with Pull Tube	Hi/Lo Motor with Pull Tube	Foot Motor with Pull Tube
5412	5492-NLA	Square 4 Pin White	7001M297	1104344-NLA	7001M324-NLA	7001M326-NLA	7001M325-NLA

## Single Drive (Careguard)

Bed Model	Foot Spring	Plug Type	Pendant (Rectangular Plug)	Head Motor with Pull Tube	Foot Motor with Pull Tube
CGBED	8050DK-NLA	Molex Plug	54000M678-NLA	7001M267-NLA	7001M263-NLA

## Electronic Conversion Kits

Conversion Type	U.S.	Canadian
Manual to Semi Electric See Note(s) A	1117169-NLA	1121066-NLA See Note(s) A
Manual to Full Electric See Note(s) B	1117170-NLA	1121068-NLA See Note(s) B
Semi to Full Electric See Note(s) C	1117171	1121067 See Note(s) C

NOTE: A - Includes: Head & Foot Motor, Junction/Control Box, Pendant, 3 Cables, Cable Lock, Motor Guard, Mounting Brackets, Motor Mounts, Cable Ties, Hardware, and Instruction Sheet.  
 B - Includes: Head, Foot, & Hi-Lo Motors, Junction/Control Box, Pendant, 4 Cables, Cable Lock, Motor Guard, Mounting Brackets, Cable Ties, Hardware, and Instruction Sheet.  
 C - Includes: Hi-Lo Motor, Junction/Control Box, Pendant, Cable, Cable Lock, Cable Ties, Hardware, and Instruction Sheet.