

Invacare Lightweight Transport LTTR/LTTB Series



Invacare Lightweight Transport LTTR/LTTB Series

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1122106-NLA	NLA - Upholstery, Back, Nylon (17" W x 18" H)	1	1
1	A	1122108-NLA	NLA - Upholstery, Back, Nylon (19" W x 18" H)	1	1
2	B	1122107-NLA	NLA - Upholstery, Seat, Nylon (17" W x 18" D)	1	1
2	B	1122109-NLA	NLA - Upholstery, Seat (19" W x 18" D)	1	1
3	C	1122102-NLA	NLA - Armrest Pad Assembly	1	2
4	D	1122101-NLA	NLA - Hardware, Hinge Plate	1	1
5		1122103-NLA	NLA - Fork Assembly, Red	1	2
5		1125020-NLA	NLA - Fork Assembly, Blue	1	2
6	E	1122105-NLA	NLA - Kit, Wheel w/ Hardware - Front	1	2
7	E	1122104-NLA	NLA - Kit, Wheel w/ Hardware - Rear	1	2
8		1122098-NLA	NLA - Cover, Step Tube - Rear	1	2
9		1122099-NLA	NLA - Guide, Seat	1	4
10		1122093-NLA	NLA - Footrest, Blue - Right	1	1
10		1122094-NLA	NLA - Footrest, Blue - Left	1	1
10		1122095-NLA	NLA - Footrest, Red - Right	1	1
10		1122096-NLA	NLA - Footrest, Red - Left	1	1
	F	1122100-NLA	NLA - Footrest Assembly, Blue (Pair)	1	1
	F	1122097-NLA	NLA - Footrest Assembly, Red (Pair)	1	1
11		1148159-NLA	NLA - Wheel Lock Assembly (Pair)	1	1
	F	1048862	Carton	1	1
	F	40131X041	Insert, Top Carton	1	1
	F	1061325-NLA	NLA - Insert, Front Removable Arm	1	1
	F	40131X039-NLA	NLA - Insert, Carton	1	2
	F	1125095-NLA	NLA - Manual, Invacare Transport, Lightweight Transport, Bariatric Transport	1	1

NOTE: A - Includes (8) Phillips Head Upholstery Screws, (2) Hex Head Screws, (2) Upper Metal Straps and (2) Lower Metal Straps.
B - Includes (8)Phillips Head Upholstery Screws, (2) Hex Head Screws and (2) Metal Straps.
C - Includes (1) Armrest Pad and Hardware.
D - Includes 4-Hinge Pins, 4 Nuts and 4 Washers
E - Includes 1-Wheel/Tire, 1-Axle Bolt and 1-Nut.
F - Not shown.