

Invacare Transport Chair with 12" Wheels

ALR/ALB Series



Invacare Transport Chair with 12" Wheels

ALR/ALB Series

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1130764-NLA	NLA - Upholstery, Back, Nylon (17" W x 18" H)	1	1
1		1130765-NLA	NLA - Upholstery, Back, Nylon (19" W x 16"D)	1	1
2		1122102-NLA	NLA - Armrest Pad Assembly	1	1
3		1122107-NLA	NLA - Upholstery, Seat, Nylon (17" W x 18" D)	1	1
3		1122109-NLA	NLA - Upholstery, Seat (19" W x 18" D)	1	1
4		1122101-NLA	NLA - Hardware, Hinge Plate	1	1
5		1122093-NLA	NLA - Footrest, Blue - Right	1	1
5		1122094-NLA	NLA - Footrest, Blue - Left	1	1
5		1122095-NLA	NLA - Footrest, Red - Right	1	1
5		1122096-NLA	NLA - Footrest, Red - Left	1	1
6		1125020-NLA	NLA - Fork Assembly, Blue	1	1
6		1122103-NLA	NLA - Fork Assembly, Red	1	1
	A	1122097-NLA	NLA - Footrest Assembly, Red (Pair)	1	1
	A	1122100-NLA	NLA - Footrest Assembly, Blue (Pair)	1	1
7		1122105-NLA	NLA - Kit, Wheel w/ Hardware - Front	1	1
8	B,D	1130766-NLA	NLA - Kit, Wheel w/ Attaching Hardware (12")	1	1
8	C	1149802-NLA	NLA - Kit, Wheel w/ Hardware - Rear (12")	1	1
9		1122099-NLA	NLA - Guide, Seat	1	1
10	E	1130767-NLA	NLA - Wheel Lock Assembly - Right	1	1
10	E	1130768-NLA	NLA - Wheel Lock Assembly - Left	1	1
	A	1130769-NLA	NLA - Wheel Lock Assembly (Pair)	1	1
	A	1122098-NLA	NLA - Cover, Step Tube - Rear	1	1
	A	1130770-NLA	NLA - Carton	1	1
	A	1130771-NLA	NLA - Insert, Top	1	1
	A	1130772-NLA	NLA - Insert, Arm	1	1
	A	1130773-NLA	NLA - Insert	1	1
	A	1125095-NLA	NLA - Manual, Invacare Transport, Lightweight Transport, Bariatric Transport	1	1
11		1132515	Kit, Hand Grip (Pair)	1	1
NOTE: A - Not shown. B - Before 10/15/07. Wheel Assembly includes wheel, bolt, locknut and hub cap. C - After 10/14/07. Wheel Assembly includes wheel, bolt, locknut and hub cap. D - Kit Part No. 1130766 is no longer available. Order Kit Part No. 1149802 as a replacement. E - When replacing the Wheel Lock Assemblies, Part Number 1132515 must be ordered also.					