



Accessories - One Arm Drive Assemblies



Accessories - One Arm Drive Assemblies

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1513269	One Arm Drive Assembly, Pneumatic w/ Hardware (15 - 17" W)	1	1
2			Shaft, Telescopic, One Arm Drive (15 - 17" W)	1	1
3			Wheel, Active, Pneumatic, DHR (24")	1	1
4			Wheel, Passive, Pneumatic, DHR (24")	1	1
5			Bushing, Quick Release	1	2
6			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	2
7			Nut, One Arm Drive	1	2
3		1509382	Wheel, Active, Pneumatic, DHR (24")	1	1
4		1508491	Wheel, Passive, Pneumatic, DHR (24")	1	1
5		1508488	Bushing, Quick Release	1	2
1	A	1513270	One Arm Drive Assembly, Pneumatic w/ Hardware (18 - 20" W)	1	1
2			Shaft, Telescopic, One Arm Drive (18 - 20" W)	1	1
3			Wheel, Active, Pneumatic, DHR (24")	1	1
4			Wheel, Passive, Pneumatic, DHR (24")	1	1
5			Bushing, Quick Release	1	2
6			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	2
7			Nut, One Arm Drive	1	2
2		1508487	Shaft, Telescopic, One Arm Drive (18 - 20" W)	1	1
3		1509382	Wheel, Active, Pneumatic, DHR (24")	1	1
4		1508491	Wheel, Passive, Pneumatic, DHR (24")	1	1
5		1508488	Bushing, Quick Release	1	2
1	A	1513271	One Arm Drive Assembly, Solid w/ Hardware (15 - 17" W)	1	1
2			Shaft, Telescopic, One Arm Drive (15 - 17" W)	1	1
3			Wheel, Solid, Active, DHR (24")	1	1
4			Wheel, Solid, Passive (24")	1	1
5			Bushing, Quick Release	1	2
6			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	2
7			Nut, One Arm Drive	1	2
2		1511769	Shaft, Telescopic, One Arm Drive (15 - 17" W)	1	1
3		1509436	Wheel, Solid, Active, DHR (24")	1	1
5		1508488	Bushing, Quick Release	1	2
1	A	1513272	One Arm Drive Assembly, Solid w/ Hardware (18 - 20" W)	1	1
2			Shaft, Telescopic, One Arm Drive (18 - 20" W)	1	1
3			Wheel, Solid, Active, DHR (24")	1	1
4			Wheel, Solid, Passive (24")	1	1
5			Bushing, Quick Release	1	2
6			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	2
7			Nut, One Arm Drive	1	2
2		1509436	Wheel, Solid, Active, DHR (24")	1	1
4		1508490	Wheel, Solid, Passive (24")	1	1

NOTE: A - Includes items 3 - 7.