


Accessories - Crutch and Cane Carrier



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	0572-02056	Cane Holder Assembly w/ Hardware	1	1
2			Strap, Hook Loop, Cane Holder	1	1
3			Holder, Cane	1	1
4			Clamp, Cane Holder	1	2
5			Spacer, Cane Clamp Holder	1	2
6			Bolt, Hex Head (M6 x 45mm)	1	2
7			NLA - Locknut (M6 x 1)	1	2
2		0572-03190	Strap, Hook Loop, Cane Holder	1	1
NOTE: A - Includes items 2 - 7.					