






# KSS Seating Accessories for Basic Head Support



## KSS Seating Accessories for Basic Head Support

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1063837-NLA	NLA - Head Support Assembly, Adult (7" D x 8" H)	1	1
1		1063838-NLA	NLA - Head Support Assembly, Child (5" D x 8" H)	1	1
2		1063840-NLA	NLA - Posts, Horizontal Head Support, Extra Long (15" L)	1	1
2		1072464-NLA	NLA - Posts, Horizontal Head Support, Adult (7" L)	1	1
2		1072442-NLA	NLA - Posts, Horizontal Head Support, Child (5" L)	1	1
3		1063839-NLA	NLA - Posts, Vertical Head Support, Extra Long (15" L)	1	1
3		1072444-NLA	NLA - Posts, Vertical, Head Support - Regular (Length 8")	1	1
4	B	0273100	Kit, Headrest Clamp, Lever Style	1	1
5	B	1072439-NLA	NLA - Knob, Adjustment, Black (3/4" L Threads)	1	1
6		1078070	Bracket Assembly, Headrest Offset, Black	1	1
	A	1063835-NLA	NLA - Pad, Head Support - Adult (3" H x 10" W)	1	1
	A	1063836-NLA	NLA - Pad, Head Support - Child (2-1/2" H x 7-1/2" W)	1	1
	A	1084842-NLA	NLA - Cover, Outer, Head Support - Adult (3" H x 10" W)	1	1
	A	1084843-NLA	NLA - Cover, Outer, Head Support - Child (2-1/2" H x 7-1/2" W)	1	1

NOTE: A - Not Shown  
 B - Item 5 Adjustment Knob was replaced with an Adjustment Lever in Headrest Clamp Kit part # 0273100. The Lever is not available as a separate part. The knob Part # 1072439 is sold separately and is interchangeable with the Lever.