


Back Cane Hardware w/o Back Brace

Adj. Height, Zero and Seven Degree



Back Cane Hardware w/o Back Brace Adj. Height, Zero and Seven Degree

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1169045-NLA	NLA - Hardware, Adjustable Height Cane w/o Backbrace	1	1
2			Bracket, Cane, Rear	1	2
3			Grip, Hand, Foam (4")	1	2
4			Insert, Backbrace	1	2
5			Screw, Socket Button Head (M6 x 53mm)	1	2
6			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
7			Screw, Cap, Socket Head (M6 x 35mm)	1	2
8			Screw, Phillips Pan Head (M4.8 X 16mm)	1	2
9			NLA - Locknut (M6 x 1)	1	6
10			Screw, Socket Head Cap (M6 x 1 x 30)	1	2
11			Panel Fastener, Snap In	1	2
2		1423410	Bracket, Cane, Rear	1	2
3		1167666	Grip, Hand, Foam (4")	1	2

NOTE: A - Includes items 2 - 11.