


Back Cane Hardware w/ Back Brace

Fixed Height, Zero and Seven Degree



Back Cane Hardware w/ Back Brace

Fixed Height, Zero and Seven Degree

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1169080	Kit, Hardware, Fixed Canes w/ Backbrace	1	1
2			Bracket, Cane, Rear	1	2
3			Grip, Cane (22 mm)	1	2
4			Insert, Rigidizer	1	2
5			Lever, Indexable	1	2
6			Screw, Cap, Socket Head (M6 x 35mm)	1	2
7			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
8			Screw, Socket Button Head (M6 x 53mm)	1	2
9			Screw, Phillips Pan Head (M4.8 X 16mm)	1	2
10			NLA - Locknut (M6 x 1)	1	2
2		1423410	Bracket, Cane, Rear	1	2
5		1514326	Lever, Indexable	1	2