



Swivel Wheel Attachments

With 13" Leg Extension and Glide Tips

Model 6269 - 3" Wheel



*Attachments Sold in Pairs
Single Side Shown*

Swivel Wheel Attachments With 13" Leg Extension and Glide Tips

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1040655	Anti-Rattle, Grey (1-1/8")	1	4
2		1040671	Tube, Leg Extension (13")	1	2
3	D	6280	Package, Tip Glide, Grey (1-1/8")	24	2
3	G	1041531-NLA	NLA - Tip, Glide, Grey (1-1/8")	1	2
3		1130781	Kit, Glide Tips	10	1
4	A,E	1039469-NLA	NLA - Caster Assembly, Swivel (3")	1	2
4	F	1123541-NLA	NLA - Caster Assembly, Swivel (3")	1	2
	B,C	21243X000-NLA	NLA - Screw, Phillips Pan Head (#6-20 x 3/8")	1	2
	B,C	1040984-NLA	NLA - Label - Right	1	1
	B,C	1040985-NLA	NLA - Label - Left	1	1
	B,C	1039493-NLA	NLA - Adapter, Bushing - Right	1	1
	B,C	1039494-NLA	NLA - Adapter, Bushing - Left	1	1

NOTE: A - Due to style changes, when replacing 3" Wheels "Before 11/5/00" order in pairs.
 B - Not Shown.
 C - Before 8/99.
 D - Includes Instruction Sheet.
 E - Item no longer available, refer to Part #1123541.
 F - Replaces Part #1039469.
 G - Individual item no longer available, refer to Kit Part #1130781.