


Front Shroud - Three Wheeler



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1517331	Reflector, Front Shroud	--	2
2		1517882	Shroud, Three Wheel Front, Silver	--	1
2		1521655	Shroud, Three Wheel Front, Blue	--	1
3		1518016	Turn Signal - Left	--	1
4		1518017	Turn Signal - Right	--	1
5		1518012	Moulding, Front Shroud	--	1
6		1527473	Head Light, Complete	--	1