


Front Frame - Four Wheeler



Front Frame - Four Wheeler

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1521732	Shroud, Front, Silver	--	1
1		1521733	Shroud, Front, Blue	--	1
1		1537792	Shroud, Front, Comet HD	--	1
2	A	1521755	Pad, Leg Shield	--	1
3	B	1521756	Pad, Front - Right/Left	--	1
4		1505476	Push Rivet, Pad	--	8
5		1505479	Tree Rivet	--	12
6		1521758	Band, Battery, Complete	--	1
7	C	1517317	Controller, Rhino (110A)	--	1
7	D	1553399	Controller, Rhino (120A)	--	1
7	C	1521759	Controller, Rhino (160A)	--	1
7	D	1553400	Controller, Rhino (160A HD)	--	1
8		1505487	Kit, Bearing, Scooter Frame	--	1
9		1521760	Front	--	1
10		1532911	Wire Harness	--	1

NOTE: A - Includes item 4.
 B - Includes item 5.
 C - For scooters manufactured prior to June 2011.
 D - For scooters manufactured after June 2011.