


# Rear Shroud



## Rear Shroud

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1521837	Rear Shroud, Silver	--	1
1		1521838	Rear Shroud, Blue	--	1
1		1530002	Shroud, Rear, Complete, Red	--	1
1		1537793	Shroud, Rear, Comet HD	--	1
2		1521862	Fender - Left	--	1
3		1521839	Fender - Right	--	1
4		1521683	Cover, Bumper, Rear	--	1
5		1518034	Turn Signal, Rear - Left	--	1
6		1518033	Turn Signal, Rear - Right	--	1
7		1521881	Reflector, Rear, E Mark	--	2
8		1518035	Neon Light, Complete	--	1
9		1521882	Cover, Rear - Left	--	1
10		1521883	Cover, Rear - Right	--	1
11		1522030	Bulb, Turn Signal (24V/10 W)	--	1
12		1522034	Bulb, Rear Light (24V/5W)	--	1