


Front Wheels - Four Wheeler



Front Wheels - Four Wheeler

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1518077	Rod Assembly, Front	--	1
2		1521994	Track Rod, Short	--	1
3		1521995	Shock Absorber, Front	--	2
4		1521996	Axle, Front - Left	--	1
5		1521997	Axle, Front - Right	--	1
6		1521999	Wheel, Front w/ Tube (13"), Complete Assembly - Grey	--	2
		1586120	Wheel, Front w/ Tube (13"), Complete Assembly - Black	--	2
		1527714	Tube (13")	--	4
		1527750	Tire, Front w/ Tube (13") - Grey	--	2
		1586118	Tire, Front w/ Tube (13") - Black	--	2
		1543159	Tire, Solid PU (13") - Grey	--	4
		1586119	Tire, Solid PU (13") - Black	--	4
7		1522000	Kit, Steering Bearing	--	2
		1527751	Kit, Tire Bearing	--	4
		1580174	Kit, Rim Front	--	2
		1509666	Kit, Mounting, Front	--	2