


# Rear Wheels and Propulsion



## Rear Wheels and Propulsion

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1522013	Transaxle (15 km/h)	--	1
1		1533710	Transaxle (12 km/h), Comet HD	--	1
2		1553429	Motor (15 km/h)	--	1
2		1533850	Motor, Comet HD	--	1
3		1505672	Lever, Motor	--	1
4		1522021	Wheel, Anti-Tipper	--	2
5		1522020	Wheel, Rear w/ Tube (13"), Complete Assembly - Grey	--	2
		1586121	Wheel, Rear w/ Tube (13"), Complete Assembly - Black	--	2
		1527714	Tube (13")	--	2
		1527750	Tire, Front w/ Tube (13") - Grey	--	2
		1586118	Tire, Front w/ Tube (13") - Black	--	2
		1543159	Tire, Solid PU (13") - Grey	--	2
		1586119	Tire, Solid PU (13") - Black	--	2
6		1533853	Brushes, Comet HD	--	4
6		1527845	Brushes, Comet	--	4
7		1533842	Guide Motor, Comet HD	--	1
8		1557459	Kit, Mounting Rear	--	4
		1505670	Washer, Plane	--	2
		1505669	Key (5 x 30)	--	2
		1527749	Kit, Rim Rear	--	2