


Front Wheel and Steering - Four Wheeler



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1518196	Spring, Front Wheel	--	2
2		1505487	Bearing, Scooter Frame	--	1
3		1518077	Rod Assembly, Front	--	1
4		1518237	Axle, Front - Left	--	1
5		1518238	Axle, Front - Right	--	1
6		1518078	Rod Assembly, Rear	--	1
7		1518333	Wheel Assembly, Front	--	2
8		1518076	Swing Assembly	--	1