


Brake Attachment (Rear Wheel)

Model 6264 - 10" Leg Extension
Model 6265 - 13" Leg Extension



*Attachments Sold in Pairs
Single Side Shown*

Brake Attachment (Rear Wheel)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1040655	Anti-Rattle, Grey (1-1/8")	1	2
2	A	1044713-NLA	NLA - Tube, Leg Extension (10")	1	2
2	B	1044929-NLA	NLA - Tube, Leg Extension (13")	1	2
3		1044944	Spring, Compression (31/64 x 1 x 3/64")	1	2
4		1044942-NLA	NLA - Housing, Spring	1	2
5		1040656	Tip, Rubber, Dark Grey (1-1/8")	1	2
6		1044943	Fork	1	2
7		7930301-NLA	NLA - Package, Button Head Cap Screw (1/4-20 x 1-3/4")	25	2
8		1039465	Bushing (1/4 x 3/8 x 1/13")	1	2
9	C	1039466-NLA	NLA - Wheel (3" x 1")	1	2
9		1130782-NLA	NLA - Package, Wheel (3")	10	1
10		7930364-NLA	NLA - Package, Locknut (1/4-20)	25	2

NOTE: A - To be used on Model 6264.
 B - To be used on Model 6265.
 C - Individual item no longer available, refer to Kit Part #1130782.