



# Walking Beam, Lower Pivot Link and Mounting Hardware

## TDX SP/SP-GT (After 3/15/11) (Canada After 6/30/11)

| Item Number | Note  | Part Number | Description  | Quantity Per |          |
|-------------|-------|-------------|--|--------------|----------|
|             |       |             |  | Package      | Assembly |
| 1           | A,C,G | 8881171997  | Kit, Walking Beam Assembly, Color, SSD, w/ Yellow Decal, Color - Right | 1            | 1        |
|             |       |             | NLA - Plug Button, Round, Black (1-1/2")                               | 1            | 1        |
|             |       |             | Decal, TDX SP, Yellow  | 1            | 1        |
|             |       |             | Screw, Socket Button Head Cap w/ Patch (1/2-13 x 7/8")                 | 1            | 1        |
| 1           | A,C,G | 8881171998  | Kit, Walking Beam Assembly, Color, SSD, w/ Yellow Decal, - Left        | 1            | 1        |
|             |       |             | NLA - Plug Button, Round, Black (1-1/2")                               | 1            | 1        |
|             |       |             | Decal, TDX SP, Yellow  | 1            | 1        |
|             |       |             | Screw, Socket Button Head Cap w/ Patch (1/2-13 x 7/8")                 | 1            | 1        |
| 1           | B,C,G | 8881171999  | Kit, Walking Beam Assembly, Color, SSD, w/ Black Decal, - Right        | 1            | 1        |
|             |       |             | NLA - Plug Button, Round, Black (1-1/2")                               | 1            | 1        |
|             |       |             | Decal, TDX SP, Black   | 1            | 1        |
|             |       |             | Screw, Socket Button Head Cap w/ Patch (1/2-13 x 7/8")                 | 1            | 1        |
| 1           | B,C,G | 8881172000  | Kit, Walking Beam Assembly, Color, SSD, w/ Black Decal, - Left         | 1            | 1        |
|             |       |             | NLA - Plug Button, Round, Black (1-1/2")                               | 1            | 1        |
|             |       |             | Decal, TDX SP, Black   | 1            | 1        |
|             |       |             | Screw, Socket Button Head Cap w/ Patch (1/2-13 x 7/8")                 | 1            | 1        |
| 3           | A     | 1144688     | Decal, TDX SP, Yellow  | 1            | 1        |
| 3           | B     | 1121689     | Decal, TDX SP, Black   | 1            | 1        |
| 5           | C,E,F | 8881142330  | Kit, Link Assembly, Lower Pivot, Color - Right                         | 1            | 1        |
| 4           |       |             | Screw, Socket Button Head Cap w/ Patch (1/2-13 x 7/8")                 | 1            | 1        |
| 6           |       |             | Screw, Socket Flat Head w/ Patch (1/2-13 x 1-1/4")                     | 1            | 1        |
| 5           | C,E,F | 8881142331  | Kit, Link Assembly, Lower Pivot, Color - Left                          | 1            | 1        |
| 4           |       |             | Screw, Socket Button Head Cap w/ Patch (1/2-13 x 7/8")                 | 1            | 1        |
| 6           |       |             | Screw, Socket Flat Head w/ Patch (1/2-13 x 1-1/4")                     | 1            | 1        |
| 7           | D     | 1145055     | Kit, Lower Pivot Link/4-Pole Motor Mounting Hardware                   | 1            | 1        |
| 8           |       |             | Locknut (1/2-20)   | 1            | 1        |
| 9           |       |             | Plate, Washer  | 1            | 1        |
| 10          |       |             | Spacer, Coved (3/4 x 1-1/8 x 1-1/32")                                  | 1            | 1        |

NOTE: A - The TDX SP yellow decal is used for all frame colors except 30P Sunny Yellow, 60P Silver Metallic and 71P Silver Vein  
B - The TDX SP black decal is used for 30P Sunny Yellow, 60P Silver Metallic and 71P Silver Vein frame colors only  
C - Specify frame color when ordering. See page titled "Frame Finishes, Touch-Up Paints and Upholstery Options".  
D - Includes items 8-10. Order two hardware kits when a total removal of the Washer Plate from the Motor Assembly "and" Lower Pivot Link Assembly is required.  
E - Lower Pivot Link for TDX SP includes the Lower Pivot Link Assembly, items 4, 6 and contains a 2 ea. press fitted plastic Flange Bearing (not shown). The Flange Bearings are not sold separately.  
F - TDX SP and TDX SC/Spree use the same Lower Pivot Link part #. When a new Lower Pivot Link is ordered, a new Screw, Socket Flat Head w/ Patch (1/2-13 x 1-1/4") will be included. The screw is used for TDX SC/Spree models only.  
G - Includes the Walking Beam Assembly and items 2-4