

Invacare Bariatric Transport Chair TRHD22FR (After 9/21/09)



Invacare Bariatric Transport Chair TRHD22FR (After 9/21/09)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1149923	NLA - Upholstery Kit, Back w/Hardware	1	1
2		1168462	Kit, Hand Grip	1	1
3		1168455-NLA	NLA - Kit, Arm Pad w/Screws (Pair)	2	1
4		1168456-NLA	NLA - Arm Assembly, Black - Right (Includes Arm, Clothing Guard, Arm Pad and Hardware)	1	1
4		1168459	Arm Assembly, Black - Left (Includes Arm, Clothing Guard, Arm Pad and Hardware)	1	1
5		1149924-NLA	NLA - Upholstery Kit, Seat w/Hardware	1	1
6		1168457	Front Rigging w/Hardware, Silver Vein - Right	1	1
6		1168458-NLA	NLA - Front Rigging w/Hardware, Silver Vein - Left	1	1
7		1168454-NLA	NLA - Kit, Wheel, Front (Includes Fork and Front Wheel Assembly, (2) Bearings, Locknut and Head Tube Cap	1	2
8		1141761	Package, Cap, End, Black, Head Tube	10	1
9		1168453-NLA	NLA - Wheel Assembly w/ Hub Cap - Rear	1	2
10		1168460-NLA	NLA - Wheel Lock Assembly - Right (Includes Complete Wheel Lock Assembly and Mounting Hardware)	1	1
10		1168461-NLA	NLA - Wheel Lock Assembly - Left (Includes Complete Wheel Lock Assembly and Mounting Hardware)	1	1
11		1149223-NLA	NLA - Hand Brake w/Cable Assembly (Pair)	2	1
	A,B	1149926-NLA	NLA - Package, Spacers	25	1
	A,B	1149927-NLA	NLA - Package, Flanged Bearings	25	1

NOTE: A - Not shown
B - Used for Rear Wheel. Bearings and Spacers are No Longer Available. A new Rear Wheel Assembly item #9 must be ordered. Includes the Hub Cap, Flanged Bearings and Spacer also.