


# Sure Glide Brakes

**Model 6274 - 10" Leg Extension**

**Model 6263 - 13" Leg Extension**



*Attachments Sold in Pairs  
Single Side Shown*

## Sure Glide Brakes

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1040655	Anti-Rattle, Grey (1-1/8")	1	2
2	A	1040847-NLA	NLA - Sure Stop Sub-Assembly (10")	1	2
2	B	1044713-NLA	NLA - Tube, Leg Extension (10")	1	2
2	B	1044929-NLA	NLA - Tube, Leg Extension (13")	1	2
3		1044942-NLA	NLA - Housing, Spring	1	2
4	A	1021734-NLA	NLA - Spring	1	2
4	B	1044944	Spring, Compression (31/64 x 1 x 3/64")	1	2
5		1044943	Fork	1	2
6		1053179	Glide	1	2
7		1040656	Tip, Rubber, Dark Grey (1-1/8")	1	2
	A,C	1021733-NLA	NLA - Tube, Inside Glide	1	2
	A,C	1040840-NLA	NLA - Plug Button, Grey (7/8")	1	2
	A,C	1041532-NLA	NLA - Tip, Glide, Grey (3/4")	1	2
NOTE: A - Before 7/95. B - After 6/95. C - Not Shown.					