


Controller, G-Trac, Programmer, 2-Pole, MK6i - Power Tiger



Controller, G-Trac, Programmer, 2-Pole, MK6i - Power Tiger

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,B	1160739-NLA	NLA - Kit, G-Trac Module Harness w/ Mounting Hardware	1	1
2			Screw, Socket Button Head (1/4-20 x 1-1/2")	1	2
3			Washer (1/4 x 5/8 x 3/32")	1	1
4			Locknut (1/4-20)	1	2
5			Locknut (#8-32)	1	2
6			NLA - Spacer, Coved (1/4 x 7/8 x 7/16")	1	1
7			NLA - Bracket, Controller Mounting, Black	1	2
8			Screw, Phillips Pan Head (#8-32 x 1")	1	2
9			Washer, Flat (3/16 x 5/8 x 1/16")	1	2
10			Harness, Module, Gyro, G-Trac	1	1
11			Screw (M6 x 1 x 20mm)	1	1
12			Nut (M6 x 1mm)	1	1
13			Bracket, Mounting, G-Trac, Black	1	1
14		1156627	Controller, G-Trac, 2-Pole, MK6i	1	1
15	C	1139985	Programmer, MK6i Remote	1	1
15	C	1142215	Kit, MK6i Remote Programmer w/ SD Professional Card & Card Holder (1813M6)	1	1
16		1139997	Cover, SD Card Slot	1	1

NOTE: A - Includes items 2-13.
 B - Item #10 part #1156626 plugs into the Controller and the other connector plugs into the Accessory Harness for Seating.
 C - The MK6i Programmer Manual (1141471) is not included and is strongly recommended the guide be ordered or downloaded from Literature Fulfillment.