


Composite Wheels - Rear w/ Hublocks, w/o Handrim

Wheel with Hub Components, Tire, Bearing
and Insert Pre-Assembled
(Hub Components Not Sold Separately)



20"
Rear Wheels

22"
Rear Wheels

24"
Rear Wheels

SCM_081A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 Rear Wheels with Hublocks	1	1
2	A	1164035	Kit, Hublock Hardware, Wheel (20" - 24")	1	1
3			Ring, Gear, Hublock	1	1
4			Screw, Button Head Cap w/ Patch (#10-24 x 1/2)	1	3

NOTE: A - Includes items 3 and 4.

Composite Wheels - Rear w/ Hublocks, w/o Handrim

#1 Rear Wheels with Hublocks

Wheel Style	Wheel Size	Tire Type	Wheel w/o Hdrn.	Tire	Tube/ Insert	Bearing
20COM	20 x 1-3/8"	Pneumatic	1164064 See Note(s) A	1017519	1017520	40117X022
20COM	20 x 1-3/8"	Pneumatic w/ Flat Free	1164065 See Note(s) A	1017519	1023322	40117X022
20COM	20 x 1-1/4"	Urethane	1164066 See Note(s) A	1037774	N/A	40117X022
22COM	22 x 1-3/8"	Pneumatic	1164067 See Note(s) A	1018544	1018545	40117X022
22COM	22 x 1-3/8"	Pneumatic w/ Flat Free	1164068 See Note(s) A	1018544	1023323	40117X022
22COM	22 x 1-1/4"	Urethane	1164069 See Note(s) A	1028250	N/A	40117X022
24COM7	24 x 1-3/8"	Pneumatic	1164070 See Note(s) A	40176X060	40176X061	40117X022
24COM7	24 x 1-3/8"	Pneumatic w/ Flat Free	1164071 See Note(s) A	40176X060	1010906	40117X022
24COM7	24 x 1-1/4"	Urethane	1164072 See Note(s) A	1026904	N/A	40117X022
24COM7	24 x 1"	KIK	1164073 See Note(s) A	1062852	N/A	40117X022

NOTE: A - Wheels come assembled with Hub Assembly, Retaining Nut, and Bearing.