


# Composite w/ Hublock and Handrim (16" and 18" Wheels)



**16COM - Rear Wheel**

Wheel w/ Handrim	Tire Type	Handrim	Handrim Hardware	Tire	Tube/ Insert	Bearing
1166050 See Note(s) A	Pneumatic	1059369	1052961	1052099	1052100	40117X022
1166049 See Note(s) A	Pneumatic w/ Flat Free	1059369	1052961	1052099	1052165	40117X022

NOTE: A - Wheel comes complete with Hub Assembly, Gear Ring, Retaining Nut, Bearing, Spring Pins, and Screws.

**18COM - Rear Wheel**

Wheel w/ Handrim	Hublock Hardware	Tire Type	Handrim	Handrim Hardware	Tire	Tube/ Insert
1104156	1164042	Pneumatic w/ Flat Free	1059369	1052961	1071441	1071140