


# Pivot/Slide Tube w/ Footplate Assembly - Composite - 1651



# Pivot/Slide Tube w/ Footplate Assembly - Composite - 1651

| Item Number | Note | Part Number | Description   | Quantity Per |          |
|-------------|------|-------------|---|--------------|----------|
|             |      |             |   | Package      | Assembly |
| 1           |      | 1026099     | Spacer, Coved (1/4 x 11/16 x 1/2")                    | 1            | 1        |
| 2           |      | 1045267     | Screw, Hex Head w/ Patch (1/4-20 x 3/4")              | 1            | 1        |
| 3           | B    | 1020534     | Insert, Threaded                                      | 1            | 1        |
| 4           |      | 1025457     | Spacer, Coved (1/4 x 3/4 x 11/32")                    | 1            | 2        |
| 5           |      | 1027097     | Plug Button, Black (3/4")                             | 1            | 1        |
| 6           |      | 60285X006   | Screw, Socket Head (1/4-20 x 2-1/4")                  | 1            | 1        |
| 7           |      | 1013519     | Washer, Nylon (5/16 x 11/16 x 1/16")                  | 1            | 1        |
| 8           |      | 1026879     | Footplate, Composite, Small, Black (6" W x 6" L)      | 1            | 1        |
| 8           |      | 1026880     | Footplate, Composite, Medium, Black (7" W x 6" L)     | 1            | 1        |
| 8           |      | 1026881     | Footplate, Composite, Large, Black (8" W x 6" L)      | 1            | 1        |
| 9           |      | 1120383     | Tube, Pivot/Slide, Curved End, Black - Short (8-1/4") | 1            | 1        |
| 9           |      | 1120382     | Tube, Pivot/Slide, Curved End, Black - Long (9-7/8")  | 1            | 1        |
| 9           |      | 1120386     | Tube, Pivot/Slide, Curved, 2" Longer, Black (11-7/8") | 1            | 1        |
| 9           |      | 1120384-NLA | NLA - Tube, Pivot/Slide, Straight - Alum. (9-7/8" L)  | 1            | 1        |
| 9           |      | 1120385     | Tube, Pivot/Slide, Straight - Black (11-7/8" L)       | 1            | 1        |
| 9           |      | 1061967     | Tube, Pivot/Slide, Straight, Black (5-1/2")           | 1            | 1        |
| 10          | A    | 1026703     | Kit, Pivot/Slide Tube Mounting Hardware (Pair)        | 1            | 1        |

NOTE: A - Item 10 includes Items 1, 2 and 3  
B - As of 8/15/03, the Threaded Insert was revised from 3 to 5 Holes. Threaded Insert is Backward Compatible and does not require ordering in Pairs.