


## Caster Forks and Hardware (After 6/29/11)



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1037420	Fork and Stem Assembly, Black (6 x 1" Caster)	1	1
1		1056640	Fork and Stem Assembly, Black (8" & 6 x 2" Caster)	1	1
1		1174562-NLA	NLA - Fork and Stem Assembly, Black (8" Caster)	1	1
2		1085298	Bearing (1-1/4" O.D. x 1/2" I.D.)	1	1
3		1001801	Bearing, Flange (1/2")	1	1
4		1021508	Washer, Nylon (5/8 x 1 x 3/32")	1	1
5		1001979	Locknut (1/2-20)	1	1
6		1001843	Cap, Head Tube, Black	1	1
6		1131189	Package, Head Tube Cap	10	1