



Caster Mounting Hardware (After 6/29/11)



Caster Mounting Hardware (After 6/29/11)

Wheel Style	Wheel Size	Tire Type	Item 1 Caster Hardware	Item 2 Axle	Item 3 Locknut	Item 4 Axle Bearing
61DC	6 x 1"	Urethane	1026926-NLA	1026908	1025730	1026160
1259	6 x 2"	Semi-Pneu.	1042489 See Note(s) A	1025726	1025734	1025724
81DC	8 x 1"	Urethane	1161851-NLA	N/A	N/A	N/A
815DC	8 x 1-1/2"	Urethane	1161851-NLA	N/A	N/A	N/A

NOTE: A - Caster Hardware #1042489 accommodates 1 side only.

NOTE: Item 1 Caster Hardware includes Axle, Locknut, and Axle Bearings.