


Caster Wheels (After 6/29/11)

Standard Caster Components

(Some Components may not be available as Stand Alone Service Parts)



Caster Wheels

Caster Model	Wheel Size	Tire Type	Item 1 Caster Wheel Assembly	Item 2 Wheel	Item 3 Tire	Item 4 Bearings	Item 5 Bearing Spacer
51DC	5"	Urethane	1169260	N/A	N/A	N/A	N/A
61DC	6"	Urethane	1164314	N/A	N/A	N/A	N/A
1259	6 x 2"	Semi- Pneumatic Urethane	1072014	N/A	N/A	1006049	1020592
81DC	8 x 1"	Urethane	1164315	N/A	N/A	N/A	N/A
815DC	8 x 1-1/2"	Urethane	1164316	N/A	N/A	N/A	N/A

NOTE: Caster Assembly includes Tire, Wheel, Bearings, and Bearing Spacer