


# Caster Wheels & Mounting Hardware (After 6/29/11)



Caster Hardware

Wheel Style	Wheel Size	Tire Type	#1 Axle	#2 Axle Nut	#3 Axle Spacer
Composite	3"/5"	Urethane	1026908	7930330-NLA	1026102
Composite	6"	Urethane	1025764	7930330-NLA	1026102
Composite	6"	Semi-Pneumatic	1025764	7930330-NLA	1026102

# Caster Wheels & Mounting Hardware (After 6/29/11)

## Caster Wheels

Wheel Style	Wheel Size	Tire Type	Bearing	Bearing Spacer	Caster Assembly Single #	Tire	Tube/ Insert
					See Note(s) A		
AT900	3"	Urethane	40117X026	1026172	1032190-NLA	N/A	N/A
41DC	4"	Urethane	N/A	N/A	1164317-NLA	N/A	N/A
51DC	5"	Urethane	N/A	N/A	1169260	N/A	N/A
6SKYWAY	6"	Urethane	40117X026	1026172	1026956	N/A	N/A
61DC	6 x 1"	Urethane	N/A	N/A	1164314	N/A	N/A
615DC	6 x 1-1/2"	Urethane	N/A	N/A	1164313	N/A	N/A

NOTE: A - Casters include Tire, Tube (if applicable), and Bearings.

6" Pneumatic Wheels with Flat-Free Inserts are No Longer Available (NLA). 6" Pneumatic Wheels are suggested as a replacement. It is recommended that they are ordered in pairs to match.