


Casters - Front (After 6/29/11)

Standard Caster Components

(Some Components may not be available as Stand Alone Service Parts)



Wheel Style	Wheel Size	Tire Type	Item 1 Caster Wheel Assembly	Item 2 Wheel	Item 3 Tire	Item 4 Bearing	Item 5 Bearing Spacer
AT900	3"	Urethane	1032190-NLA	N/A	N/A	40117X026	1026172
51DC	5 x 1"	Urethane	1169260	N/A	N/A	N/A	N/A
61DC	6 x 1"	Urethane	1164314	N/A	N/A	N/A	N/A
615DC	6 x 1-1/2"	Urethane	1164313	N/A	N/A	N/A	N/A
1259 (After 8/14/03)	6 x 2" Deep Profile	Semi Pneumatic	1072014	N/A	N/A	1006049	1020592
81DC	8 x 1"	Urethane	1164315	N/A	N/A	N/A	N/A
815DC	8 x 1-1/2"	Urethane	1164316	N/A	N/A	N/A	N/A

NOTE: Caster Wheel Assembly includes Wheel, Tire, Bearings (2), and Bearing Spacer.