








Casters and Hardware



Casters and Hardware

| Item Number | Note | Part Number | Description | Quantity Per | |
|-------------|------|-------------|--|--------------|----------|
| | | | | Package | Assembly |
| 1 | | 1159772 | Kit, Front Wheel Hardware | 1 | 1 |
| | | | Spacer (5/16 x 7/16) | 1 | 2 |
| | | | Locknut (5/16-24) | 1 | 1 |
| | | | Screw, Hex Head (5/16-24 x 2-1/4") | 1 | 1 |
| | | | See Chart Caster Styles | 1 | 1 |
| 2 | | 1128070 | Wheel Assembly, Soft Roll (4" x 1") | 1 | 1 |
| 2 | | 1128071 | Wheel Assembly, Soft Roll (5" x 1-1/2") | 1 | 1 |
| 2 | | 1157282-NLA | NLA - Wheel Assembly, Soft Roll (3" x 1-1/2") | 1 | 1 |
| 2 | | 1169260 | Wheel Assembly, Composite, Urethane (5" x 1") | 1 | 1 |
| 3 | | 1032190-NLA | NLA - Wheel Assembly (3" x 1") | 1 | 1 |
| 4 | | 1060910 | Wheel Assembly, Composite Urethane, Light Up (4" x 1") | 1 | 1 |
| 5 | | 1164317-NLA | NLA - Wheel Assembly, Composite, Urethane (4 x 1") | 1 | 1 |
| 5 | | 1164318-NLA | NLA - Wheel Assembly, Composite, Urethane (4 x 1-1/2") | 1 | 1 |
| 6 | | 1114142 | Wheel Assembly, Semi-Pneumatic, Black (4" x 1") | 1 | 1 |
| 7 | | 1164311 | Wheel Assembly, Composite, Urethane (5" x 1") | 1 | 1 |

Caster Styles

| Caster Style | Size | Tire Type | Part Number |
|--------------|--------|--------------------|-------------|
| 3SOFR | 3" | Aluminum Soft Roll | 1157282-NLA |
| CRFF3 | 3" | Micro | 1032190-NLA |
| 41DC | 4 x 1" | Urethane | 1164317-NLA |
| CRFA654 | 4" | Aluminum Soft Roll | 1128070 |
| CRF844 | 4" | Lighted | 1060910 |
| 51DC | 5 x 1" | Urethane | 1169260 |
| CRFA655 | 5" | Aluminum Soft Roll | 1128071 |