


Arms - T Style Adjustable - Myon Active (REDL, REFL)



Arms - T Style Adjustable - Myon Active (REDL, REFL)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1449987	T-Arm Assembly, Central w/ Hardware, Desk Length - Right	1	1
2			Pad, Arm, Desk Length	1	1
3			Arm Pad Mounting Assembly - Left	1	1
4			Tube, Mounting, Adjustable	1	1
5			NLA - Support, Armrest	1	1
6			Guard, Clothing	1	1
7			Bracket, Clothing Guard	1	2
8			Screw, Socket Head (M12 x 12.5)	1	1
9			Bushing, Armrest	1	1
10			Nut, Square, Armrest	1	1
11			Screw, Countersunk (M4.8)	1	2
12			Screw (M4.2 - 16 8.8)	1	4
13			Screw, Socket Head (M6 x 50)	1	2
14			NLA - Locknut (M6 x 1)	1	2
1		1449988	T-Arm Assembly, Central w/ Hardware, Desk Length - Left	1	1
2			Pad, Arm, Desk Length	1	1
3			Arm Pad Mounting Assembly - Right	1	1
4			Tube, Mounting, Adjustable	1	1
5			NLA - Support, Armrest	1	1
6			Guard, Clothing	1	1
7			Bracket, Clothing Guard	1	2
8			Screw, Socket Head (M12 x 12.5)	1	1
9			Bushing, Armrest	1	1
10			Nut, Square, Armrest	1	1
11			Screw, Countersunk (M4.8)	1	2
12			Screw (M4.2 - 16 8.8)	1	4
13			Screw, Socket Head (M6 x 50)	1	2
14			NLA - Locknut (M6 x 1)	1	2
1		1450015	T-Arm Assembly, Central w/ Hardware, Full Length - Right	1	1
2			Pad, Arm, Full Length	1	1
3			Arm Pad Mounting Assembly - Left	1	1
4			Tube, Mounting, Adjustable	1	1
5			NLA - Support, Armrest	1	1
6			Guard, Clothing	1	1
7			Bracket, Clothing Guard	1	2
8			Screw, Socket Head (M12 x 12.5)	1	1
9			Bushing, Armrest	1	1
10			Nut, Square, Armrest	1	1
11			Screw, Countersunk (M4.8)	1	2
12			Screw (M4.2 - 16 8.8)	1	4
13			Screw, Socket Head (M6 x 50)	1	2
14			NLA - Locknut (M6 x 1)	1	2
1		1450016	T-Arm Assembly, Central w/ Hardware, Full Length - Left	1	1
2			Pad, Arm, Full Length	1	1

Arms - T Style Adjustable - Myon Active (REDL, REFL)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
3			Tube, Mounting, Adjustable	1	1
4			Arm Pad Mounting Assembly - Right	1	1
5			Screw, Socket Head (M6 x 50)	1	2
6			NLA - Support, Armrest	1	1
7			Guard, Clothing	1	1
8			Bracket, Clothing Guard	1	2
9			Screw, Socket Head (M12 x 12.5)	1	1
10			Bushing, Armrest	1	1
11			Nut, Square, Armrest	1	1
12			Screw, Countersunk (M4.8)	1	2
13			Screw (M4.2 - 16 8.8)	1	4
14			Screw, Socket Head (M6 x 50)	1	2