


# Adductors Mounting Hardware, LNX/LNX-L Front Riggings

For replacement Adductors, see section titled, "Accessories-Adductors"



# Adductors Mounting Hardware, LNX/LNX-L Front Riggings

For replacement Adductors, see section titled, "Accessories-Adductors"

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1162434	Hardware, Mounting, Adductor, MPCTMT/MPCTMT-L Front Riggings	1	1
2			Plug Button, Black (1")	1	2
NOTE: A - Includes item 2. Adductor Mounting Hardware used for MPCTMT/MPCTMT-L and LNX/LNX-L Power Center Mount Articulating Elevating Legrest to plug the end of the Seat Rails only.					