

Wheels - Rear w/o Handrims - Spoke (After 9/1/11)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	A	1030936	Plug, Valve Stem	1	1
	A	See Chart	Wheels without Handrims - Spoke	1	1

NOTE: A - Not Shown

Wheels - Rear w/o Handrims - Lightweight Performance (ITESPK)

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube Insert	Bearing	Hub Spacer
22ITESPK	22"	Pneumatic	1176799-NLA	1018544	1018545	N/A	N/A
22ITESPK	22"	Pneumatic w/ Flat Free	1176800-NLA	1018544	1023323	N/A	N/A
24ITESPK	24"	Pneumatic	1176802	1026138	1032240	N/A	N/A
26ITESPK	26"	Pneumatic	1176803-NLA	N/A	N/A	N/A	N/A

Wheels - Rear w/o Handrims - Standard Performance (SP)

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube Insert	Bearing	Hub Spacer
24SP	24"	Continental	1116705-NLA See Note(s) A	1026139-NLA	N/A	40117X021	1015756
24SP	24"	Pneumatic	1115660 See Note(s) A	40176X060	40176X061	40117X021	1015756
24SP	24"	Pneumatic w/ Flat Free	1116707 See Note(s) B, A	40176X060	1010906	40117X021	1015756
24SP	24"	KIK	1114112 See Note(s) A	1062852	N/A	40117X021	1015756
26SP	26"	Continental	1116721-NLA See Note(s) A	1026143-NLA	N/A	40117X021	1015756
26SP	26"	Primo	1116723-NLA See Note(s) A	1029470	1029471	40117X021	1015756

NOTE: A - Stem Hole Plugs (Part Number 1030936) are used with Flat Free and Urethane Tires.
 B - As of (11/08/02) these wheels now have a Black finish. For chairs built Before (11/08/02), it is recommended when one wheel is requested that the wheels be replaced in pairs.

Wheels - Rear w/o Handrims - Spoke (After 9/1/11)

Wheels - Rear w/o Handrims - High Performance (HPS)

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube Insert	Bearing	Hub Spacer
24HPS	24"	Continental	1116713-NLA See Note(s) A	1026139-NLA	N/A	40117X021	1015756
24HPS	24"	Pneumatic	1116709-NLA See Note(s) B, A	40176X060	40176X061	40117X021	1015756
24HPS	24"	Pneumatic (Pair)	1116710	40176X060	40176X061	40117X021	1015756
24HPS	24"	Pneumatic w/ Flat Free	1116715 See Note(s) D, A	40176X060	1010906	40117X021	1015756
24HPS	24"	Primo	1116711 See Note(s) A	1026138	1032240	40117X021	1015756
24HPS	24"	KIK	1114138 See Note(s) A	1062852	N/A	40117X021	1015756
26HPS	26"	Continental	1115653-NLA See Note(s) C	1026143-NLA	N/A	40117X021	1015756
26HPS	26"	Pneumatic	1115652-NLA See Note(s) C	1026142	1026140	40117X021	1015756
26HPS	26"	Pneumatic	1116725-NLA	N/A	N/A	N/A	N/A
26HPS	26"	Primo	1115654 See Note(s) A	1029470	1029471	40117X021	1015756

NOTE: A - As of (11/08/02) these wheels now have a Black finish. For chairs built Before (11/08/02), it is recommended when one wheel is requested that the wheels be replaced in pairs.
 B - Wheels must be replaced in pairs. See Part number (1116710) for replacement.
 C - No longer available as an individual service part. Wheel must be replaced in pairs, use 1116725 as its replacement.
 D - Stem Hole Plugs (1030936) are used with Flat Free and Urethane Tires.

Wheels - Rear w/o Handrims - Spoke (After 9/1/11)

Wheels - Rear w/o Handrims - Spinergy (SPIN)

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube Insert	Bearing	Hub Spacer
SPIN24	24"	Primo (Single #)	1112338-NLA See Note(s) B, A	1026138	1105941	40117X021	N/A
SPIN24	24"	Primo (1 Pair)	1112353-NLA See Note(s) B, A	1026138	1105941	1104350	N/A
SPIN25	25"	Primo (Single #)	1112348-NLA See Note(s) B, A	1044947	1104326	1104350	N/A
SPIN25	25"	Primo (1 Pair)	1112354-NLA See Note(s) B, A	1044947	1104326	1104350	N/A
SPIN26	26"	Primo (Single #)	1112343-NLA See Note(s) B, A	1029470	1104326	1104350	N/A
SPIN26	26"	Primo (1 Pair)	1112355-NLA See Note(s) B, A	1029470	1104326	1104350	N/A

NOTE: A - As of (7/15/02) Spinergy Wheels now have 18 spokes instead of 30. For chairs built Before (7/15/02), it is recommended when one wheel is requested that the wheels be replaced in pairs.
B - As of (11/17/08) Spinergy Wheels changed from black spoke, black hub to a black spoke, silver hub appearance. For chairs built Before (11/17/08), it is recommended when one wheel is requested that the wheels be replaced in pairs.

Wheels - Rear w/o Handrims - Miscellaneous (559TS, 700SP)

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube Insert	Bearing	Hub Spacer
559TS	25"	Primo	1116717-NLA See Note(s) A	1044947	1029471	40117X021	1015756
700SP	27"	Umma Gumma	1116727-NLA See Note(s) A	1045836-NLA	1045837-NLA	40117X021	1015756

NOTE: A - As of (11/08/02) these wheels now have a Black finish. For chairs built Before (11/08/02), it is recommended when one wheel is requested that the wheels be replaced in pairs.