

Spoke Rear Wheels (w/ Handrims) (After 9/1/11)

Wheel Style	Tire Type	Anodized Aluminum	Grey Plastic Coated	Black Plastic Coated	8 Vertical Projection	8 Oblique Projection
20ITESPK	Pneumatic	1176796	N/A	N/A	N/A	N/A
20SP	Pneumatic w/ Flat Free	1115586 See Note(s) A	1049450-NLA	1115424-NLA See Note(s) A	N/A	N/A
22ITESPK	Pneumatic	1176801	N/A	N/A	N/A	N/A
22SP	Pneumatic with Flat Free	1115590-NLA See Note(s) A	1049456-NLA	1115593 See Note(s) A	1115426-NLA See Note(s) A	N/A
24SP	Pneumatic	1115594-NLA See Note(s) A	1049460-NLA	1115599 See Note(s) A	1115427-NLA See Note(s) A	1115428-NLA See Note(s) A
24SP	Pneumatic with Flat Free	1115597 See Note(s) A	1049465-NLA	1115601 See Note(s) A	1115429-NLA See Note(s) A	1115430-NLA See Note(s) A
24SP	Primo	1115602-NLA See Note(s) A	1049470-NLA	1115603-NLA See Note(s) A	1115431-NLA See Note(s) A	1115432-NLA See Note(s) A
24SP	Continental	1115606-NLA See Note(s) A	1049475-NLA	1115607-NLA See Note(s) A	1115433-NLA See Note(s) A	1115434-NLA See Note(s) A
24SP	KIK	1115645 See Note(s) A	1091332-NLA	1115646 See Note(s) A	1115447-NLA See Note(s) A	1115448-NLA See Note(s) A

NOTE: A - As of 11/8/02 these wheels now have a Black finish. For chairs built Before 11/8/02, it is recommended when one wheel is requested that the wheels be replaced in pairs.

When ordering a Spoke Wheel with Handrim attached the standard handrim mounting process uses the outside holes of the handrim mounting tabs. To use the inside mounting holes to reduce the overall width of the chair, you must order the wheel and handrim assembled that way as a special order.