


Composite Rear Wheel w/ Hublock (16" Wheel)



Wheel Type	Wheel Size	Tire Type	Item 1 Wheel Assembly	Item 2 Tire	Item 3 Tube/Insert
16COM	16" x 1-3/4"	Pneumatic	1166046 See Note(s) A	1052099	1052100
16COM	16" x 1-3/4"	Pneumatic w/ Flat Free	1166047 See Note(s) A	1052099	1052165

NOTE: A - Wheel comes complete with Hub Assembly, Gear Ring, Retaining Nut, Bearing, Spring Pins, Tire, Insert, and Screws.