


Basket, Floor Mats and Accessories

Leo 3W, Leo 4W



Basket, Floor Mats and Accessories Leo 3W, Leo 4W

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	B	1164840	Floor Mats w/ Hardware, Leo 4 Wheel	1	1
2	C	1547589	Floor Mats w/ Hardware, Leo 3 Wheel	1	1
3	A,D	1505414-NLA	NLA - Basket	1	1
3	A,E	1148861-NLA	NLA - Basket	1	1
3	A	1577430	Front Composite Basket	1	1
4		1164841	Package, Push Rivets	25	1
5		1164837	Reflector (Pair)	2	1

NOTE: A - Additional hardware is not needed to attach basket.
 B - Includes floor mats and attaching hardware. For Leo 4 Wheel Scooter only.
 C - Includes floor mat and attaching hardware. For Leo 3 Wheel Scooter only.
 D - Part number no longer available. See part number 1148861 as a replacement.
 E - Part number no longer available. See part number 1577430 as a replacement.