


Seat & Arms

Leo 3W, Leo 4W



Seat & Arms Leo 3W, Leo 4W

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1173792-NLA	NLA - Arm Assembly - Right	1	1
1	A	1173793-NLA	NLA - Arm Assembly - Left	1	1
1	B	1506314	Arm Assembly - Right	1	1
1	B	1506313	Arm Assembly - Left	1	1
2		1164833	Seat Slide Assembly	1	1
3		1164835-NLA	NLA - Seat and Back Assembly	1	1
4		1164834-NLA	NLA - Seat Bracket Assembly	1	1
5		1165106	Post, Seat Guide w/ Hardware	1	1
6		1173789	Armrest	1	2
7	C	1168463-NLA	NLA - Cover, Back, Leo - Black Only	1	1
8	C	1168464-NLA	NLA - Cover, Seat, Leo - Black Only	1	1

NOTE: A - For Four Wheel Scooters manufactured before 08/08/2011.
 B - For Three or Four Wheel Scooters manufactured after 08/07/2011.
 C - For complete Seat Assembly, see item 3 (1164835).