


Matrix PCS Back Headrest Clamp Kit - Lever Style

(Matrix Brand name replaced InTouch as of 6/30/2011)



Matrx PCS Back Headrest Clamp Kit - Lever Style

(Matrx Brand name replaced InTouch as of 6/30/2011)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1160332	Kit, Headrest Clamp, Knob Style	1	1
2		1072439-NLA	NLA - Knob, Adjustment, Black (3/4" L Threads)	1	1
3	B	0237522-NLA	NLA - Clamp, Small Top Half	1	1
4		1025736	Locknut (1/4-20)	1	2
5	B	0237422-NLA	NLA - Clamp, Small Bottom Half	1	1
6		1009838	Screw, Socket Button (1/4-20 x 3/4")	1	2
NOTE: A - Includes items 2-6. B - Clamp is no longer available individually. See Headrest Clamp Kit item #1.					