


# Matrix PCS Lateral Support

(Matrix Brand name replaced InTouch as of 6/30/2011)



## Matrix PCS Lateral Support

(Matrix Brand name replaced InTouch as of 6/30/2011)

Item Number	Note	Part Number	Description	Quantity Per		
				Package	Assembly	
1	A	1164438	Kit, Lateral Support w/ Small Pad (Pair)	1	1	
2			NLA - Screw, Socket Head (1/4-20 x 1/2")	1	4	
3			Washer (1/4 x 1/2 x 1/16")	1	4	
4			Lateral, w/ Foam - Small	1	2	
5			Upholstery Assembly, Lateral w/ Label - Small	1	1	
5			Upholstery Assembly, Lateral - Small	1	1	
1	A	1164439	Kit, Lateral Support w/ Large Pad (Pair)	1	1	
2			NLA - Screw, Socket Head (1/4-20 x 1/2")	1	4	
3			Washer (1/4 x 1/2 x 1/16")	1	4	
4			Lateral, w/ Foam - Large	1	2	
5			Upholstery Assembly, Lateral w/ Label - Large	1	1	
5			Upholstery Assembly, Lateral - Large	1	1	
4			1164820	Lateral, w/ Foam - Large	1	2
4			1164821	Lateral, w/ Foam - Small	1	2
5	B	1165480	Upholstery Assembly, Lateral w/ Label - Small	1	2	
5	B	1165482	Upholstery Assembly, Lateral w/ Label - Large	1	2	

NOTE: A - Includes items 2-5.  
 B - Upholstery Assembly w/o a Label is not saleable.