


# Wheel Locks - Pull to Lock

Pull to Lock (WLPULL), Pull to Lock Extended (WLPULLX)



## Wheel Locks - Pull to Lock

### Pull to Lock (WLPULL), Pull to Lock Extended (WLPULLX)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1538079	Wheel Lock Assembly w/ Hardware, Pull to Lock - Right	1	1
1		1538078	Wheel Lock Assembly w/ Hardware, Pull to Lock - Left	1	1
2		1541957	Wheel Lock Assembly w/ Hardware, Pull to Lock w/ Extensions - Right	1	1
2		1541956	Wheel Lock Assembly w/ Hardware, Pull to Lock w/ Extensions - Left	1	1
3	A	1423107	Kit, Hardware, Clamp Brake (Pair)	1	1
4			Screw (M6 x 1 x 35mm)	1	4
5			Mounting Clamp, Wheel Lock	1	2

NOTE: A - Includes items 4 and 5.