

Wheels - Rear w/o Handrims - Spoke (After 9/1/11)

Wheels - Rear w/o Handrims - Lightweight Performance (ITESPK)

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube/Insert	Bearing	Hub Spacer
20ITESPK	20"	Pneumatic	1176797	1017519	1017520	N/A	N/A
20ITESPK	20"	Pneumatic w/ Flat Free	1176798	1017519	1023322	N/A	N/A
22ITESPK	22"	Pneumatic	1176799-NLA	1018544	1018545	N/A	N/A
22ITESPK	22"	Pneumatic w/ Flat Free	1176800-NLA	1018544	1023323	N/A	N/A
24ITESPK	24"	Pneumatic	8881159332 See Note(s) A	40176X060	1032240	N/A	N/A
24ITESPK	24"	Pneumatic w/ Flat Free	8881159332 See Note(s) A	40176X060	1010906	N/A	N/A
24ITESPK	24"	Primo	8881159332 See Note(s) A	1026138	1032240	N/A	N/A
24ITESPK	24"	KIK	8881159332 See Note(s) A	1062852	N/A	N/A	N/A
24ITESPK	24"	Schwalbe	8881159332 See Note(s) A	1139813	1105941	N/A	N/A
26ITESPK	26"	Pneumatic	1176803-NLA	1026142	1026140	1026141-NLA	N/A
26ITESPK	26"	Primo	N/A	1029470	1029471	N/A	N/A
26ITESPK	26"	Schwalbe	1176804	1139815	1104326	N/A	N/A

NOTE: A - Use the Service Parts Wheel Configurator when ordering a complete wheel for this style. To omit handrim in the Service Parts Wheel Configurator, order as a special.

Wheels - Rear w/o Handrims - Spoke (After 9/1/11)

Wheels - Rear w/o Handrims - High Performance (HPS)

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube/ Insert	Bearing	Hub Spacer
24HPS	24 x 1-3/8"	Pneumatic	1116709-NLA See Note(s) A	40176X060	40176X061	40117X021	1015756
24HPS	24 x 1-3/8"	Pneumatic (Pair)	1116710	40176X060	40176X061	40117X021	1015756
24HPS	24 x 1-3/8"	Pneumatic w/ Flat Free	1116715	40176X060	1010906	40117X021	1015756
24HPS	24 x 1"	Primo	1116711	1026138	1032240	40117X021	1015756
24HPS	24 x 1"	KIK	1114138	1062852	N/A	40117X021	1015756
24HPS	24 x 1"	Schwalbe	1157361-NLA	1139813	1105941	40117X021	N/A
26HPS	26 x 1-3/8"	Pneumatic	1115652-NLA	1026142	1026140	40117X021	1015756
26HPS	26 x 1-3/8"	Pneumatic	1116725-NLA	N/A	N/A	N/A	N/A
26HPS	26 x 1"	Primo	1115654	1029470	1029471	40117X021	1015756
26HPS	26 x 1"	Schwalbe	1154389-NLA	1139815	1104326	40117X021	N/A

NOTE: A - Wheels must be replaced in pairs. See Part number 1116710 for replacement.

Wheels - Rear w/o Handrims - Spoke (After 9/1/11)

Wheels - Rear w/o Handrims - Spinergy (SPIN)

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube/Insert	Bearing	Hub Spacer
SPIN24 Black Spokes	24 x 1"	Primo	1112338-NLA See Note(s) A	1026138	1105941	1104350	N/A
SPIN24 Yellow, Red, Blue or White Spokes	24 x 1"	Primo	8881159332 See Note(s) B	1026138	1105941	1104350	N/A
SPIN24 Black Spokes	24 x 1"	Schwalbe	1154229	1139813	1105941	1104350	N/A
SPIN24 Yellow, Red, Blue or White Spokes	24 x 1"	Schwalbe	8881159332 See Note(s) B	1139813	1105941	1104350	N/A
SPIN25 Black Spokes	25 x 1"	Primo	1112348-NLA	1044947	1104326	1104350	N/A
SPIN25 Yellow, Red, Blue or White Spokes	25 x 1"	Primo	8881159332 See Note(s) B	1044947	1104326	1104350	N/A
SPIN25 Black Spokes	25 x 1"	Primo Treaded Pneumatic	1130593-NLA	1130527	1104326	1104350	N/A
SPIN25 Yellow, Red, Blue or White Spokes	25 x 1"	Primo Treaded Pneumatic	8881159332 See Note(s) B	1130527	1104326	1104350	N/A
SPIN25 Black Spokes	25 X 1"	Schwalbe	1154232-NLA	1139814	1104326	1104350	N/A
SPIN25 Yellow, Red, Blue or White Spokes	25 x 1"	Schwalbe	8881159332 See Note(s) B	1139814	1104326	1104350	N/A
SPIN26 Black Spokes	26 x 1"	Primo	1112343-NLA	1029470	1104326	1104350	N/A
SPIN26 Yellow, Red, Blue or White Spokes	26 x 1"	Primo	8881159332 See Note(s) B	1029470	1104326	1104350	N/A
SPIN26 Black Spokes	26 x 1"	Schwalbe	1154235-NLA	1139815	1104326	1104350	N/A
SPIN26 Yellow, Red, Blue or White Spokes	26 x 1"	Schwalbe	8881159332 See Note(s) B	1139815	1104326	1104350	N/A

NOTE: A - As of 11/17/08 Spinergy Wheels changed from black spoke, black hub to a black spoke, silver hub appearance. For chairs built before 11/17/08, it is recommended when one wheel is requested that the wheels be replaced in pairs.
B - Use the Parts Wheel Configurator when ordering a complete wheel for this model. To omit handrim in the Parts Wheel Configurator, order as a special

Wheels - Rear w/o Handrims - Spoke (After 9/1/11)

Wheels - Rear w/o Handrims - Spinergy LX (SPLX)

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube/Insert	Bearing	Hub Spacer
SPLX24 Black, Yellow, Red, Blue or White Spokes	24 x 1"	Primo	8881159332 See Note(s) A	1026138	1105941	1104350	N/A
SPLX24 Black, Yellow, Red, Blue or White Spokes	24 x 1"	Schwalbe	8881159332 See Note(s) A	1139813	1105941	1104350	N/A
SPLX25 Black, Yellow, Red, Blue or White Spokes	25 x 1"	Primo	8881159332 See Note(s) A	1044947	1104326	1104350	N/A
SPLX25 Black, Yellow, Red, Blue or White Spokes	25 x 1"	Primo Treaded Pneumatic	8881159332 See Note(s) A	1130527	1104326	1104350	N/A
SPLX25 Black, Yellow, Red, Blue or White Spokes	25 x 1"	Schwalbe	8881159332 See Note(s) A	1139814	1104326	1104350	N/A

NOTE: A - Use the Parts Wheel Configurator when ordering a complete wheel for this model. To omit handrim in the Parts Wheel Configurator, order as a special

Wheels - Rear w/o Handrims - Sun Fusion (FN)

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube/Insert	Bearing	Hub Spacer
FN1624	24 x 1"	Primo	8881159332 See Note(s) A	1026138	1105941	40117X021	N/A
FN1624	24 x 1"	Pneumatic	8881159332 See Note(s) A	40176X060	40176X061	40117X021	N/A
FN1624	24 x 1"	KIK	8881159332 See Note(s) A	1062852	N/A	40117X021	N/A
FN1624	24 x 1"	Schwalbe	8881159332 See Note(s) A	1139813	1105941	40117X021	N/A
FN1625	25 x 1"	Primo	8881159332 See Note(s) A	1044947	1104326	40117X021	N/A
FN1625	25 x 1"	Primo Treaded Pneumatic	8881159332 See Note(s) A	1130527	1104326	40117X021	N/A
FN1625	25 x 1"	KIK	8881159332 See Note(s) A	1144435	N/A	40117X021	N/A
FN1625	25 x 1"	Schwalbe	8881159332 See Note(s) A	1139814	1104326	40117X021	N/A
FN1626	26 x 1"	Primo	8881159332 See Note(s) A	1029470	1104326	40117X021	N/A

Wheels - Rear w/o Handrims - Spoke (After 9/1/11)

Wheels - Rear w/o Handrims - Sun Fusion (FN)

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube/Insert	Bearing	Hub Spacer
FN1626	26 x 1"	Pneumatic	8881159332 See Note(s) A	1026142	1026140	40117X021	N/A
FN1626	26 x 1"	Schwalbe	8881159332 See Note(s) A	1139815	1104326	40117X021	N/A
NOTE: A - Use the Parts Wheel Configurator when ordering a complete wheel for this model. To omit handrim in the Parts Wheel Configurator, order as a special							

Wheels - Rear w/o Handrims - Spoke (After 9/1/11)

Wheels - Rear w/o Handrims - 559 Turbo Spoke (559TS)

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube/Insert	Bearing	Hub Spacer
559TS	25 x 1"	Primo	1116717-NLA	1044947	1029471	40117X021	1015756
599TS	25 x 1"	Primo Treaded Pneumatic	1130592-NLA	1130527	1104326	40117X021	1015756