


# Seat Frame Mounting Hardware & Interface Mtg. Plates (CGR) For TDX SP/SP-GT/SR w/ 3-22NF/GP24 Batteries After 7/13/11



## Seat Frame Mounting Hardware & Interface Mtg. Plates (CGR) For TDX SP/SP-GT/SR w/ 3-22NF/GP24 Batteries After 7/13/11

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1143803-NLA	NLA - Kit, Seat Mounting Plate Hardware (Recline Only)	1	1
2			Screw, Hex Head (1/4-20 x 7/8")	1	8
3			Locknut (1/4-20)	1	8
4	B	1147142	Kit, Formula CGT/CGR/CGTR Interface Asbly. Mtg. Hdw. (For TDX SP w/ GP 24 Batteries Less Vent Option)	1	1
5			Screw, Hex Head w/ Patch (5/16-18 x 3/4")	1	4
6			Washer, Flat (21/64 x 53/64 x 7/64")	1	4
4	C	1147146-NLA	NLA - Kit, Formula CG Tilt, Recline & Tilt/Recline Interface Asbly. Mtg. Hdw. (For TDX SP w/ GP 24 Batteries and Vent Option)	1	1
5			NLA - Screw, Socket Head (5/16-18 x 3/4")	1	4
6			Washer, Flat (21/64 x 53/64 x 7/64")	1	4
7	C	1147130	Kit, Interface Mounting Plate Hardware	1	1
8			Plate, Nut	1	2
9			Screw, Hex Head SEMS w/ Patch (5/16-18 x 3/4")	1	8
10	D	1173239	Plate, Low Interface Mountng, Black - Right	1	1
10	D	1173240	Plate, Low Interface Mountng, Black - Left	1	1
10	E	1173241	Plate, Medium Interface Mountng, Black - Right	1	1
10	E	1173242	Plate, Medium Interface Mountng, Black - Left	1	1
10	F	1173243	Plate, High Interface Mountng, Black - Right	1	1
10	F	1173244	Plate, High Interface Mountng, Black - Left	1	1
NOTE: A - Includes items 2 and 3 B - Includes items 5 and 6 C - Includes items 8 and 9 D - Low Seat-To-Floor = 17.25-18.25" adjustment range. E - Medium Seat-To-Floor = 18-19" adjustment range. F - High Seat-To-Floor = 19-20" adjustment range.					