

Front Riggings Reference Chart

Pivot/ Slide tube and footplates are interchangeable with the support assemblies listed above them.

1651=Composite footplate with 9-7/8" curved pivot/slide tube.
AT5543=Adjustable angle footplate with 9-7/8" straight pivot/slide tube.
1450=Deluxe aluminum footplate with "L" style pivot/slide tube.
1350=Extra large footplate with "L" style pivot/slide tube.
NOTE: Calf pads with attaching hardware not included on legrest.
NOTE: Pivot/ Slide tube attaching hardware not included.
NOTE: Pivot/Slide tube and footplates do not include heel loops. See footrest and legrest accessories.

Support Assemblies

Rigging Type	AT5544	60HD	70NHD	PW93	PW93E	60 MFX
Right	8881045291 See Note(s) A	8881053158 See Note(s) A	8881053149 See Note(s) A	8881070776 See Note(s) A	8881070768-NL A See Note(s) A	8881026998 See Note(s) A
Left	8881045292 See Note(s) A	8881053157 See Note(s) A	8881053148 See Note(s) A	8881070777 See Note(s) A	8881070769-NL A See Note(s) A	8881026999 See Note(s) A
NOTE: A - Specify frame color when ordering.						

Pivot/ Slide Tube with Footplates

Footplate Type	1651 13-14" W	1651 15-16" W	1651 17-24" W	1350 17-24" W	AT5543 12-13" W	AT5543 14-15" W	AT5543 16-24" W
Right	1057859	1057869	1057879	1026395	1060442	1060444	1062756
Left	1057860	1057870	1057880	1026396	1060443	1060445	1062757
Right 2" Longer	1057861	1057871	1057881	1086515	N/A	1062758	1062760
Left 2" Longer	1057862-NLA	1057872	1057882	1086517	N/A	1062759	1062761
Right 4" Longer	1058419-NLA	1058421	1058423	1086516	N/A	N/A	N/A
Left 4" Longer	1058420-NLA	1058422	1058424	1086518	N/A	N/A	N/A

Front Riggings Reference Chart

Support Assembly

AT5544

Right	Left
8881045291	8881045292
See Note(s) A	See Note(s) A
NOTE: A - Specify frame color when ordering	

Pivot/ Slide Tube and Footplates

1450 16"-24" Wide

Right	Left
1028723-NLA	1028724-NLA

STD. = Aluminum footplate with "L" style pivot/slide tube.
1450 = Deluxe aluminum footplate with "L" style pivot/slide tube.
1350 = Extra large footplate with "L" style pivot/slide tube.
NOTE: Calf pads with attaching hardware not included on legrest.
NOTE: Pivot/Slide tube attaching hardware not included.
NOTE: Pivot/Slide tube and footplates do not include heel loops. See footrest and legrest accessories.

Support Assembly

P904A

Right	Left
8881058646-NLA	8881058647-NLA
See Note(s) A	See Note(s) A
NOTE: A - Specify frame color when ordering.	

Front Riggings Reference Chart

Pivot/ Slide Tube and Footplates

P904A

Footplate Type	Standard 14-15" W	Standard 16-24" W	1450 16-24" W	1350 17-24" W
Right	1058612-NLA	1058604	1028723-NLA	1058830
Left	1058613-NLA	1058605	1028724-NLA	1058831
Right 2" Longer	N/A	1086529-NLA	N/A	1086541-NLA
Left 2" Longer	N/A	1086531-NLA	N/A	1086551-NLA
Right 4" Longer	N/A	1086530-NLA	N/A	1086542-NLA
Left 4" Longer	N/A	1086532-NLA	N/A	1086552-NLA

Support Assembly

PAL4A/PAS4A (Chrome Only)

Right	Left
1058482-NLA	1058483-NLA

Pivot/ Slide Tube and Footplates

PALA4/PAS4A

Footplate Type	Standard 14-15" W	Standard 16-24" W	1450 16-24" W	1350 17-24" W
Right	1023905-NLA	1023907	1023924-NLA	1035456
Left	1023906-NLA	1023908	1023925-NLA	1035457

Front Riggings Reference Chart

1651=Composite footplate with 9-7/8" curved pivot/slide tube.
AT5543 = Adjustable angle footplates with 9-7/8" straight pivot/slide tube.
Std. = Aluminum footplate with "L" style pivot/slide tube.
1450 = Deluxe aluminum footplate with "L" style pivot/style tube.
1350 = Extra large footplate with "L" style pivot/slide tube.
5043 = Flip up footplates without pivot/slide tubes.
NOTE: Pivot/Slide tube attaching hardware not included.
NOTE: Pivot/Slide tube and footplates do not include heel loops. See footrest and legrest accessories.

Support Assembly

P93M

Right	Left
8881037836-NLA	8881037837-NLA
See Note(s) A	See Note(s) A
NOTE: A - Specify frame color when ordering	

Pivot/ Slide Tube and Footplates

P93M

Footplate Type	Standard 14-15" W	Standard 16-24" W	1450 16-24" W	1350 17-24" W
Right	1058573-NLA	1057969-NLA	1028892-NLA See Note(s) A	1031924-NLA
Left	1058574-NLA	1057970-NLA	1028893-NLA See Note(s) A	1031925-NLA
Right 2" Longer	N/A	N/A	N/A	1086519-NLA
Left 2" Longer	N/A	N/A	N/A	1086523-NLA
Right 4" Longer	N/A	N/A	N/A	1086520-NLA
Left 4" Longer	N/A	N/A	N/A	1086524-NLA
NOTE: A - No longer available as a service part individually.				

Front Riggings Reference Chart

70 MFX and 90 Swing Supports

Supports	70 MFX	90 Swing
Right	8881047480 See Note(s) A	8881062596 See Note(s) A
Left	8881047481 See Note(s) A	8881062597 See Note(s) A
NOTE: A - Specify frame color when ordering.		

Footplate Assemblies

70 MFX and 90 Swing Supports

Footplate Type	AT5043 12-16" W	AT5543 12-13" W	AT5543 14-15" W	AT5543 17-24" W
Right	1029434-NLA	1060436-NLA	1060438	1062570
Left	1029435-NLA	1060437	1060439	1062571

Support Assembly

PW93ST/ 70STTAPER

Rigging Type	PW93ST	70STTAPER
Right	8881070784-NLA See Note(s) A	8881026996 See Note(s) A
Left	8881070785-NLA See Note(s) A	8881026997 See Note(s) A
NOTE: A - Specify frame color when ordering.		

Front Riggings Reference Chart

Pivot/Slide Tube and Footplates

PW93ST / 70STTAPER

Footplate Type	1651 15-17" W	1651 18-19" W	1651 20-24" W
Right	1057859	1057869	1057879
Left	1057860	1057870	1057880
Right 2" Longer	1057861	1057871	1057881
Left 2" Longer	1057862-NLA	1057872	1057882