


# Wheels - Drive Wheels for Non-Free Wheel Hubs

Gearless/Brushless motors were made available on-chair as an option for Storm Arrow, Ranger X and Torque on 11/2/98. On 7/1/99, Gearless/Brushless motors became standard for Storm Arrow only.



**Item #s 1, 2, and 3**

Wheel Size	Item #1 Complete Wheel	Item #1 Complete Wheel w/ Flat Free	Item #2 Tire	Item #2 Tire w/ Flat Free	Item #3 Tube
14 x 4" 4 Pole Motors	1042476	1042477	1042474	1042478	1042475
14 x 4" GB Motors  See Note(s) A	1079305-NLA	1079304	1042474	1042478	1042475
14 x 3" 4 Pole and GB Motors	1039653	1041449-NLA	1039652	1101754	1039651
NOTE: A - For drive wheel rims and attaching hardware for the Gearless/Brushless 14 x 4" Wheel, see "Drive Wheels and Attaching Hardware for Free Wheel Hub Mounts".					

## Wheels - Drive Wheels for Non-Free Wheel Hubs

Gearless/Brushless motors were made available on-chair as an option for Storm Arrow, Ranger X and Torque on 11/2/98. On 7/1/99, Gearless/Brushless motors became standard for Storm Arrow only.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1-3		See Chart	#s 1, 2, and 3 Drive Wheels for Non-Free Wheel Hubs	1	1
4		1042453-NLA	NLA - Kit, 6" Rim w/ Hardware (14" x 4" Wheel)	1	1
4		1090651-NLA	NLA - Kit, 8" Rim w/ Hardware (14" x 3" Wheel)	1	1
	A	1008030-NLA	NLA - Screw, Hex Head (3/8-16 x 3/4")	1	6
	A	1005519	NLA - Lockwasher (3/8 x 11/16 x 3/32")	1	6
	A	18058X000-NLA	NLA - Nut, Hex (3/8-16)	1	6
	B	1024443-NLA	NLA - Screw, Hex Head (5/16-24 x 5/8")	1	6
	B	1024442-NLA	NLA - Washer, Lock (5/16 x 9/16 x 1/16")	1	6
	B	1024426-NLA	NLA - Nut, Hex (5/16-24)	1	6
	C,D	7930326-NLA	NLA - Package, Screw, Hex Head (1/4-20 x 3/4")	25	8
	C,D	7930333	Locknut (1/4-20)	25	8
	D	40176X010-NLA	NLA - Pump, Tire	1	1
NOTE: A - Included with Part Number 1042453 before 3/4/98. B - Included with Part Number 1042453 after 3/3/98. C - Included with Part Number 1090651. D - Not Shown.					