



Arms - T Style Adjustable - Myon (TAHD, TAHF)



Arms - T Style Adjustable - Myon (TAHD, TAHF)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1428596-0601	Armrest, T-Arm, Adjustable Height, Desk - Right	1	1
3			Arm Pad, Desk Length	1	1
5			Tube, Armrest - Right	1	1
6			Guard, Clothing - Right/Left	1	1
7			NLA - Support, Armrest	1	1
8			Screw, Socket, Button Head (M5 x 35mm)	1	2
10			Plug Button, Black (3/4")	1	2
11			Screw, Socket Head (M6 x 50)	1	2
12			Screw, Socket Button Head (M6 x 1 x 12mm)	1	3
13			NLA - Locknut (M6 x 1)	1	5
16			Button, Spring	1	1
16		0027-10003	Button, Spring	1	1
1	A	1428597-0601	Armrest, T-Arm, Adjustable Height, Desk - Left	1	1
3			Arm Pad, Desk Length	1	1
5			Tube, Armrest, Black - Left	1	1
6			Guard, Clothing - Right/Left	1	1
7			NLA - Support, Armrest	1	1
8			Screw, Socket, Button Head (M5 x 35mm)	1	2
10			Plug Button, Black (3/4")	1	2
11			Screw, Socket Head (M6 x 50)	1	2
12			Screw, Socket Button Head (M6 x 1 x 12mm)	1	3
13			NLA - Locknut (M6 x 1)	1	5
16			Button, Spring	1	1
16		0027-10003	Button, Spring	1	1
2	B	1428598-0601	Armrest, T-Arm, Adjustable Height, Full - Right	1	1
4			Arm Pad, Full Length	1	1
5			Tube, Armrest - Right	1	1
6			Guard, Clothing - Right/Left	1	1
7			NLA - Support, Armrest	1	1
8			Screw, Socket, Button Head (M5 x 35mm)	1	2
9			Button, Plug (22mm)	1	1
10			Plug Button, Black (3/4")	1	2
11			Screw, Socket Head (M6 x 50)	1	2
12			Screw, Socket Button Head (M6 x 1 x 12mm)	1	3
13			NLA - Locknut (M6 x 1)	1	5
14			Screw, Socket Head, Torx (K50 x 30mm)	1	1
15			Tube, Support, Armpad, Full Length	1	1
16			Button, Spring	1	1
16		0027-10003	Button, Spring	1	1
2	B	1428599-0601	Armrest, T-Arm, Adjustable Height, Full - Left	1	1
4			Arm Pad, Full Length	1	1
5			Tube, Armrest, Black - Left	1	1
6			Guard, Clothing - Right/Left	1	1
7			NLA - Support, Armrest	1	1
8			Screw, Socket, Button Head (M5 x 35mm)	1	2
9			Button, Plug (22mm)	1	1
10			Plug Button, Black (3/4")	1	2

Arms - T Style Adjustable - Myon (TAHD, TAHF)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
11			Screw, Socket Head (M6 x 50)	1	2
12			Screw, Socket Button Head (M6 x 1 x 12mm)	1	3
13			NLA - Locknut (M6 x 1)	1	5
14			Screw, Socket Head, Torx (K50 x 30mm)	1	1
15			Tube, Support, Armpad, Full Length	1	1
16			Button, Spring	1	1
16		0027-10003	Button, Spring	1	1
9		1151389	Package, Plug Button, Black (3/4")	10	1
10		1151390-NLA	NLA - Package, Plug Button, Black (22mm)	10	1
NOTE: A - Includes items 3, 5 - 8, 10 - 13, and 16. B - Includes items 4 - 16.					