


Seat and Arms - L-3X

NOTE: A full seat assembly is not available at this time, please order seat components separately.



Seat and Arms - L-3X

NOTE: A full seat assembly is not available at this time, please order seat components separately.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1168468	Arm Assembly (Pair)	2	1
			Arm Assembly - Right	1	1
			Arm Assembly - Left	1	1
2	C	1148847-NLA	NLA - Cover, Seat	1	1
2	D	1168466	Cover, Seat	1	1
3	C	1148848-NLA	NLA - Cover, Backrest	1	1
3	D	1168465	Cover, Backrest	1	1
4		1148845	Kit, Seat Slide	1	1
5		1148846-NLA	NLA - Seat Bracket Assembly	1	1
5	D	1169832	Seat Bracket Assembly	1	1
6		1148850-NLA	NLA - Hardware, Seat Post	1	1
			Includes: (1) ea. Nylon Nut, (1) ea. Bolt, & (2) ea. Washers.		
7		1148852	Liner, Shroud/Guide Seat	1	1
8		1148849-NLA	NLA - Kit, Plugs and Knobs	1	1
			Includes: (2) ea. Screws & (2) ea. Plugs.		
9		1144812	Spring, Seat Lever	1	1

NOTE: C - Before 9/7/10
D - After 9/6/10